

**FREE TRIAL**[Home](#) / [FAQs](#) / Question and Answer

How big are the FreeStyle Libre Sensors? How deeply is it inserted?

The FreeStyle Libre 2 Sensor is 5mm in height and 35mm in diameter.

The FreeStyle Libre 3 Sensor is even smaller at just 2.9mm in height and 21mm in diameter.

The portion of both Sensors that gets inserted under the skin is less than 0.4 mm wide and it's inserted only about 5mm under the skin. In a study conducted by Abbott Diabetes Care, users agreed that inserting the Sensor was no more painful than a typical finger prick glucose test*.

*Data on file. In a study conducted by Abbott Diabetes Care, 93.4% of patients surveyed (n=30) did not feel discomfort under the skin while wearing the Sensor.

[RETURN TO THE FAQ PAGE](#)[DISCOVER](#)[PRODUCTS](#)[GETTING STARTED](#)[MYFREESTYLE](#)[SUPPORT](#)

[Privacy Notice](#)[Terms and Conditions of Use](#)[Terms of Sale](#)[Cookie Policy](#)[Cookie Preferences](#)

The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners. No use of any Abbott trademark, trade name, or trade dress in this site may be made without the prior written authorisation of Abbott Laboratories, except to identify the product or services of the company.

This website and the information contained herein is intended for use by residents of the United Kingdom. The product images are for illustrative purposes only.

© 2024 Abbott Laboratories Limited. All rights reserved.

Registered Number: 329102 England. Registered Office: Abbott House, Vanwall Business Park, Vanwall Road, Maidenhead, Berkshire, SL6 4XE