



US010941438B2

(12) **United States Patent**
Wang et al.

(10) **Patent No.: US 10,941,438 B2**
(45) **Date of Patent: Mar. 9, 2021**

(54) **METHODS AND KITS FOR HIGHLY
MULTIPLEX SINGLE PRIMER EXTENSION**

FOREIGN PATENT DOCUMENTS

(71) Applicant: **QIAGEN SCIENCES, LLC,**
Germantown, MD (US)

WO 98/01562 A1 1/1998
WO 2007/123742 A2 11/2007
WO 2013/175815 A1 11/2013

(72) Inventors: **Yexun Wang,** Ellicott City, MD (US);
Quan Peng, Clarksburg, MD (US)

OTHER PUBLICATIONS

(73) Assignee: **QIAGEN Sciences, LLC,** Germantown,
MD (US)

(*) Notice: Subject to any disclaimer, the term of this
patent is extended or adjusted under 35
U.S.C. 154(b) by 367 days.

"Bisulfite sequencing," Wikipedia, archived Nov. 2, 2017, URL=
[https://en.wikipedia.org/w/index.php?title=Bisulfite_sequencing](https://en.wikipedia.org/w/index.php?title=Bisulfite_sequencing&oldid=808370619)
&oldid=808370619, download date Jul. 11, 2018, 12 pages.
Altschuler, "PCR Troubleshooting," URL=[http://www.highveld.com/](http://www.highveld.com/pcr/pcr-troubleshooting.html)
pcr/pcr-troubleshooting.html, download date Jul. 11, 2018, 5 pages.
Biswas et al., "Identification and Characterization of a Thermostable
MutS Homolog from *Thermus aquaticus*," *The Journal of Biological*
Chemistry 271(9):5040-5048, 1996. (10 pages).
Brotherton et al., "Single primer extension (SPEX) amplification to
accurately genotype highly damaged DNA templates," *Forensic*
Science International: Genetics Supplement Series 1(1):19-21, 2008.
Eisen, "A phylogenomic study of the MutS family of proteins,"
Nucleic Acids Research 26(18):4291-4300, 1998.
Fukui et al., "Thermostable Mismatch-Recognizing Protein MutS
Suppresses Nonspecific Amplification during Polymerase Chain
Reaction (PCR)," *Int. J. Mol. Sci.* 14:6436-6453, 2013.
Mitani et al., "Rapid SNP diagnostics using asymmetric isothermal
amplification and a new mismatch-suppression technology," *Nature*
Methods 4(3):257-262, 2007.
Myllykangas et al., "Efficient targeted resequencing of human
germline and cancer genomes by oligonucleotide-selective sequenc-
ing," *Nature Biotechnology* 29(11):1024-1027, 2011. (26 pages).
Peng et al., "Reducing amplification artifacts in high multiplex
amplicon sequencing by using molecular barcodes," *BMC Genom-
ics* 16:589, 2015. (12 pages).
Primrose et al., *Principles of Gene Manipulation and Genomic*, 7th
ed., Blackwell Publishing, Hoboken, New Jersey, USA, 2005, Chap.
8, "Changing genes: site-directed mutagenesis and protein engi-
neering," pp. 141-156.
Rizzi et al., "Ancient DNA studies: new perspectives on old
samples," *Genetics Selection Evolution* 44:21, 2012. (19 pages).
Stanislawska-Sachadyn et al., "MutS as a tool for mutation detec-
tion," *Acta Biochimica Polonica* 52(3):575-583, 2005.
Watanabe et al., "Use of a competitive probe in assay design for
genotyping of the UGT1A1*28 microsatellite polymorphism by the
smart amplification process," *BioTechniques* 43(4):479-484, 2007.

(21) Appl. No.: **15/767,312**

(22) PCT Filed: **Oct. 14, 2016**

(86) PCT No.: **PCT/US2016/057073**

§ 371 (c)(1),

(2) Date: **Apr. 10, 2018**

(87) PCT Pub. No.: **WO2017/066592**

PCT Pub. Date: **Apr. 20, 2017**

(65) **Prior Publication Data**

US 2019/0062822 A1 Feb. 28, 2019

Related U.S. Application Data

(60) Provisional application No. 62/242,766, filed on Oct.
16, 2015.

(51) **Int. Cl.**

C12Q 1/68 (2018.01)

C12N 9/12 (2006.01)

C12Q 1/6848 (2018.01)

C12Q 1/6844 (2018.01)

C12Q 1/6869 (2018.01)

C12Q 1/686 (2018.01)

(52) **U.S. Cl.**

CPC **C12Q 1/6848** (2013.01); **C12N 9/1252**
(2013.01); **C12Q 1/686** (2013.01); **C12Q**
1/6844 (2013.01); **C12Q 1/6869** (2013.01)

(58) **Field of Classification Search**

None

See application file for complete search history.

(56) **References Cited**

U.S. PATENT DOCUMENTS

5,877,280 A 3/1999 Wetmur
6,150,112 A 11/2000 Weissman et al.
6,287,825 B1 9/2001 Weissman et al.
2008/0044921 A1 2/2008 Iwaki et al.

Primary Examiner — James Martinell

(74) Attorney, Agent, or Firm — Seed IP Law Group LLP

(57) **ABSTRACT**

The present disclosure provides methods and kits for highly
multiplex single primer extensions using a MutS protein and
Mg2+ at a concentration higher than that in a typical PCR
reaction. Also disclosed is the use of such methods and kits
in next generation sequencing.



25 Claims, 4 Drawing Sheets


Specification includes a Sequence Listing.


EXHIBIT #

Metzker 2212

12/03/2025

1. DNA random fragmentation**2. Library construction with NGSCounter****3. Single primer extension****4. Universal PCR library****5. Library quantification****6. Sequencing****FIG. 1**

**FIG. 2**

**FIG. 3**




FIG. 4

1

METHODS AND KITS FOR HIGHLY MULTIPLEX SINGLE PRIMER EXTENSION

CROSS REFERENCE TO RELATED APPLICATION

This application claims the benefit under 35 U.S.C. § 119(e) to U.S. Provisional Application No. 62/242,766, filed Oct. 16, 2015, which application is hereby incorporated by reference in its entirety.

STATEMENT REGARDING SEQUENCE LISTING

The Sequence Listing associated with this application is provided in text format in lieu of a paper copy, and is hereby incorporated by reference into the specification. The name of the text file containing the Sequence Listing is 830109_409USPC_SEQUENCE_LISTING.txt. The text file is 37.8 KB, was created on Mar. 15, 2018, and is being submitted electronically via EFS-Web.

BACKGROUND

Technical Field

The present disclosure relates to highly multiplex single primer extensions with improved efficiency and specificity and their use in next generation sequencing (NGS).

Description of the Related Art

Targeted sequencing is one of the major applications for NGS. It is a cost efficient way to sequence the most relevant portion of the genome or transcriptomes.

Highly multiplex PCR is one of the popular methods to prepare samples for targeted sequencing. Although it is a simple yet powerful tool to enrich a large gene region, it has limited ability to detect large structural variants such as gene fusion or splicing variants because the two primers per amplicon design cannot accommodate large or unknown structural changes between the two predefined primer sites. In addition, highly multiplex PCR is difficult to be scaled up to cover very large regions such as whole transcriptome or exome in a single reactions tube. The many pairs of primers will introduce unwanted artifacts when multiplex scale increases dramatically.

A recent alternative approach to enriching target loci is single primer extensions. The single primer design allows the detection of large and unknown structural variants, and improves scalability for covering very larger panels. However, because target specificity is provided by only one primer, the overall enrichment specificity is generally much lower than PCR enrichment using two target-specific primers.

One way to improve single primer specificity is to conduct the single primer annealing at more stringent conditions (e.g., at high temperature and high salt condition) and only perform polymerase extension after primer binding has reached equilibrium. However, high stringent hybridization conditions typically reduce DNA polymerase activity. In addition, performing polymerase extension after primer binding is inconvenient.

Another challenge in enriching a very large region is the limitation of each primer concentration. Because of the physical and chemical constraints, when more different primers are added in one reaction (to cover more regions),

2

the concentration of each primer inevitably has to be lowered in proportion. For example, each primer may need to be at 2 nM when 20,000 primers are used to enrich 500 genes, while each primer may be at 20 nM when 2000 primers are used to enrich 50 genes. In other cases, the concentration is limited by how the primer is synthesized. For example, oligonucleotides prepared by massive parallel synthesis on microchip have much lower yields than by traditional column-based synthesis approach. In both scenarios, lower primer concentrations will lead to lower template annealing efficiency under standard conditions.

One way to improve annealing efficiency is to prolong the annealing time to hours, which is not very desirable. While it is possible to change the annealing chemistry to accelerate proper primer-template binding, doing so increases non-specific primer binding and extension by DNA polymerase, leading to overall low specificity.

SUMMARY

In one aspect, the present disclosure provides a method for performing a highly multiplex single primer extension reaction, comprising:

a. extending at least about 1,000 different primers in a single primer extension reaction using a plurality of target nucleic acids as templates in the presence of (1) a MutS protein, and (2) Mg^{2+} at a minimum concentration of about 6 mM to generate extension products.

In certain embodiments, at least about 4,000 additional different primers, such as at least 9,000 additional different primers are extended in the single primer extension reaction.

In certain embodiments, the single primer extension reaction is performed in the absence of a primer capable of specifically annealing to a portion of the extension products.

In certain embodiments, most of the target nucleic acids comprise a common sequence, and wherein the reaction is performed in the presence of a boosting primer that comprises the common sequence or a portion thereof that is at least about 10 nucleotides in length.

In certain embodiments, each of the primers other than the boosting primer is present in the single primer extension reaction at a concentration of no more than about 20 nM.

In certain embodiments, each of the primers other than the boosting primer is present in the single primer extension reaction at a concentration of about 1 nM to about 10 nM.

In certain embodiments, Mg^{2+} is present in the single primer extension reaction at a concentration of about 6 mM to about 20 mM.

In certain embodiments, the plurality of target nucleic acids are genomic DNA or amplification products of genomic DNA.

In certain other embodiments, the plurality of target nucleic acids are cDNA or amplification products of cDNA.

In certain embodiments, the plurality of target nucleic acids are bisulfite treated genomic DNA or its amplified products.

In certain embodiments, the MutS protein is a *Thermus aquaticus* MutS protein.

In certain embodiments, the primer annealing and extension specificity of the single primer extension reaction is at least about 80%.

In certain embodiments, the average primer annealing and extension efficiency of the single primer extension reaction is at least about 10%.

In certain embodiments, the method provided herein further comprises:

b. optionally amplifying the extension product(s) to generate amplified extension product(s), and

c. sequencing the extension product(s) or the amplified extension products.

In certain embodiments, the target nucleic acids are partially double-stranded, and comprise a double-stranded target nucleic acid sequence and a single-stranded common sequence at the 5' terminus of each strand of the target nucleic acid sequence.

In certain embodiments, prior to step a. the method disclosed herein comprises:

(i) fragmenting DNA molecules to generate double-stranded DNA fragments, and

(ii) ligating the single-stranded common sequence to the 5' terminus of each strand of the double-stranded DNA fragment to provide partially double-stranded target nucleic acids.

In certain embodiments, each primer other than the boosting primer if present comprises at its 5' terminus a universal sequence that is incapable of specifically annealing to a region of a target nucleic acid.

In a related aspect, the present disclosure provides a method for sequencing nucleic acids, comprising:

a. providing a plurality of partially double-stranded target nucleic acids each of which comprises a double-stranded target nucleic acid sequence, and a single-stranded common sequence at the 5' terminus of each strand of the target nucleic acid sequence, wherein the single-stranded common sequence is the same among the plurality of the target nucleic acids;

b. providing a plurality of primers capable of specifically annealing to portions of the target nucleic acid sequences;

c. extending the primers using the target nucleic acid sequences as templates in the presence of (1) a MutS protein and (2) Mg^{2+} at a minimum concentration of about 6 mM to obtain extension products; and

d. sequencing the extension products or their amplification products.

In certain embodiments, extending the primers in step c. is performed in the presence of a boosting primer comprising the common sequence or a portion thereof that is at least about 10 nucleotides in length.

In certain embodiments, step a. comprises:

a1. fragmenting double-stranded DNA molecules to generate double-stranded DNA fragments, and

a2. ligating the single-stranded common sequence to the 5' terminus of each strand of the double-stranded DNA fragments to provide the plurality of partially double-stranded nucleic acids.

In certain embodiments, the method comprises: prior to step d., a further step of amplifying the extension products of step c.

In another aspect, the present application provides a kit for performing a primer extension reaction, comprising:

(1) a pool of at least about 1,000 different primers,

(2) a buffer concentrate comprising a MutS protein, dNTPs, and Mg^{2+} , wherein the concentration of Mg^{2+} in a highly multiplex single primer extension reaction mixture generated from diluting the buffer concentrate is at least about 6 mM, and

(3) a DNA polymerase.

In certain embodiments, the kit comprises a pool of at least about 10,000 different primers.

In certain embodiments, the MutS protein is a *Thermus aquaticus* MutS protein.

BRIEF DESCRIPTION OF THE DRAWINGS

FIG. 1 is a schematic representation of an embodiment of the present disclosure, showing that single primer extension is used in sequencing.

FIG. 2 is a graph showing primer annealing and extension efficiencies at various Mg^{2+} concentrations.

FIG. 3 is a graph comparing primer annealing and extension specificities between an extension reaction performed using 2 nM each primer and 12 mM Mg^{2+} and an extension reaction performed using 20 nM each primer and 4 mM Mg^{2+} .

FIG. 4 is a graph comparing primer annealing and extension specificities and uniformities (T50) between an extension reaction performed at 12 mM Mg^{2+} without MutS and an extension reaction performed at 12 mM Mg^{2+} with MutS.

DETAILED DESCRIPTION

Targeted sequence enrichment using single primer extension approach has great advantages over traditional two primer PCR amplicon based enrichment. Because one boundary of the amplified product is not constrained by one predefined primer, it has the ability to detect unexpected structural variations in the sequence, such as fusions and splicing variants. However, due to only one primer being used, single primer extension has traditionally suffered from low enrichment specificity. Another challenge in target enrichment is the low primer annealing efficiency especially when the primer concentration is low to accommodate large gene panels. Existing approaches to promoting primer binding typically inadvertently reduce binding specificity.

Faced with the challenges of low primer annealing efficiency and specificity, the present inventors have identified a solution to increasing primer annealing extension efficiency while improving primer extension specificity via a combination of Mg^{2+} at a concentration higher than in a typical PCR reaction and a MutS protein. This combination is especially beneficial when a large number of very low concentrations of primers are used in single primer extension process to enrich a large genomics region for sequencing analysis.

In the following description, any ranges provided herein include all the values in the ranges. It should also be noted that the term "or" is generally employed in its sense including "and/or" (i.e., to mean either one, both, or any combination thereof of the alternatives) unless the content dictates otherwise. Also, as used in this specification and the appended claims, the singular forms "a," "an," and "the" include plural referents unless the content dictates otherwise. The terms "include," "have," "comprise" and their variants are used synonymously and to be construed as non-limiting. The term "about" refers to $\pm 10\%$ of a reference value. For example, "about 50° C." refers to "50° C. $\pm 5^\circ$ C." (i.e., 50° C. $\pm 10\%$ of 50° C.).

I. Methods for Performing Single Extension Reactions

In one aspect, the present disclosure provides a method for performing a highly multiplex single primer extension reaction that comprises extending at least about 1,000 different primers in a single primer extension reaction using a plurality of target nucleic acids as templates in the presence of (1) a MutS protein, and (2) Mg^{2+} at a minimum concentration of about 6 mM to generate extension products.

The term "single primer extension" or "single primer extension reaction" as used in the present disclosure refers to a reaction in which a primer is extended using a single-stranded target nucleic acid or one strand of a double-

stranded target nucleic acid as a template. A single primer extension does not include extension of another primer that anneals to the complementary strand of the single-stranded target nucleic acid or the other strand of the double-stranded target nucleic acid except where a boosting primer is used as described below.

A “primer” is an oligonucleotide capable of specifically annealing to a target nucleic acid, which allows the extension of the primer using the target nucleic acid as a template. Generally, such a primer or a sufficiently long portion thereof (e.g., at least about 10, 11, 12, 13, 14, or 15 nucleotides long) is substantially or completely complementary to a region the target nucleic acid to allow specific annealing of the primer to the target nucleic acid. A primer used in a multiplex single primer extension reaction as disclosed herein for targeted nucleic acid enrichment may also be referred to as an “enrichment primer.”

An “oligonucleotide” refers to a short polymer composed of deoxyribonucleotides, ribonucleotides or combinations thereof. Oligonucleotides are generally between about 10 to about 100 nucleotides, preferably about 12 to about 60 nucleotides, in length.

The terms “complementary” and “complement” and their variants, as used herein, refer to any two nucleic acid sequences or portions thereof that form a hybridized duplex by base pairing. One nucleic acid fragment may be completely complementary to another nucleic acid fragment if all of the nucleotides in the fragment form base pairing with nucleotides in the corresponding antiparallel positions on the other nucleic acid fragment. A primer (or a portion thereof that is at least about 10 nucleotides in length) is “substantially complementary” to a target nucleic acid if at least 90% (e.g., at least 95%, at least 98%, or at least 99%) of the residues of the primer (or the portion thereof) are complementary to corresponding residues in the target nucleic acid over the length of the primer (or the portion thereof).

In certain embodiments, a primer is substantially or completely complementary to a target nucleic acid over the complete length of the primer. In certain other embodiments, a primer is not substantially complementary to a target nucleic acid over the complete length of the primer. In such embodiments, the primer may comprise a sequence at its 3' end that is substantially or completely complementary to the target nucleic acid so that the primer may specifically anneal to the target nucleic acid, and another sequence at its 5' end that is not substantially complementary to the target nucleic acid sequence.

A single primer extension reaction may include multiple cycles of annealing of a primer to a single-stranded target nucleic acid or a strand of a double-stranded target nucleic acid, extending the primer, and denaturing the extension product from its complementary template to increase the amount of extension product.

In certain embodiments, the single primer extension reaction does not comprise any primer capable of specifically annealing to the complementary strand of the single-stranded target nucleic acid or the other strand of the double-stranded target nucleic acid. In such embodiments, multiple cycles of annealing, extending, and denaturing lead to linear amplification of the extension product.

A primer is capable of specifically annealing to a template nucleic acid if the primer is substantially or completely complementary to the template nucleic acid or a portion thereof that is at least about 10 nucleotides in length.

In certain other embodiments, the single primer extension reaction comprises a boosting primer. As discussed below, including a boosting primer in a single primer extension

reaction that comprises multiple cycles of denaturing, annealing and extending will lead to exponential amplification of the extension product.

A “multiplex single primer extension reaction” is a single primer extension reaction where a plurality of primers (also referred to as “a primer set” or “a pool of primers”) other than a boosting primer as described below are extended using one or more target nucleic acids as templates.

A “highly multiplex single primer extension reaction” refers to a multiplex single primer extension reaction where at least about 1,000 different primers other than a boosting primer are extended using a plurality of target nucleic acids as templates.

In certain embodiments, a highly multiplex single primer extension reaction may contain about 1,000 to about 100,000 different primers. For example, a highly multiplex single primer extension may contain at least about 1,000, 1,500, 2,000, 2,500, 3,000, 3,500, 4,000, 4,500, or 5,000 different primers, and/or at most about 10,000, 20,000, 30,000, 40,000, 50,000, 60,000, 70,000, 80,000, 90,000, or 100,000 different primers, including any combination of the above-noted minimum and maximum numbers of different primers such as about 1,000 to about 10,000, about 1,000 to about 50,000, about 1,000 to about 100,000, about 5,000 to about 20,000, about 5,000 to about 50,000, or about 5,000 to about 100,000 different primers.

In a highly multiplex single primer extension reaction, primers other than a boosting primer (if present) may anneal to different genes, different regions of the same genes, and/or different strands of the same regions of the same genes as long as none of such primers forms a pair with another primer in the same reaction that allows exponentially amplification of a portion of a target nucleic acid delineated by the primer pair. In cases where two primers in a multiplex single primer extension reaction anneal to different portions of the same strand of a target nucleic acid, preferably, the different portions are sufficiently apart from each other so that the extension of one primer does not interfere with the extension of the other primer.

In certain embodiments, each primer other than a boosting primer (if present) in a highly multiplex single primer extension reaction comprises at its 5' terminus a universal sequence that is incapable of specifically annealing to a region of a target nucleic acid. Such a universal sequence is useful in designing additional primers for amplifying the products of the highly multiplex single primer extension reaction. Additional description of such a universal sequence is provided in Section II below.

The primers may be obtained in any methods known in the art for synthesizing oligonucleotides. In certain embodiments, the primers are synthesized on microchips.

In a multiplex single primer extension reaction, the initial concentration (i.e., the concentration at the start of the single primer extension reaction) of each primer other than a boosting primer may be about 0.01 nM to about 100 nM, preferably about 1 nM to about 40 nM, such as about 2 nM to about 20 nM.

In certain embodiments, each primer other than a boosting primer in a multiplex single primer extension reaction has an initial concentration of at least about 0.01, 0.02, 0.03, 0.04, 0.05, 0.06, 0.07, 0.08, 0.09, 0.1, 0.2, 0.3, 0.4, or 0.5 nM and/or at most about 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 30, 40, 50, 60, 70, 80, 90, or 100 nM, including any combination of the above-noted minimum and maximum concentrations.

Target nucleic acids useful as templates in single primer extension reactions disclosed herein may be any nucleic

acids of interest, including genomic DNA, cDNA, genomic DNA treated with bisulfite, and the amplification products of the above-noted nucleic acids.

Target nucleic acids may be single-stranded, double-stranded, or partially double-stranded. In certain embodiments, target nucleic acids are partially double-stranded and comprise a double-stranded target nucleic acid sequence and a single-stranded common sequence at the 5' terminus of each strand of the target nucleic acid sequence. The single-stranded common sequence is the same among the target nucleic acids, and allows the design of a boosting primer as described below to increase or maximize the yield of extension products.

Target nucleic acids may be isolated from any samples, sources, or organisms of interest that contain nucleic acids, including materials obtained from clinical, forensic, and environmental settings. Preferably, the sample is a biological sample derived from a human, animal, plant, bacteria or fungi. The term "sample" also includes processed samples including preserved, fixed and/or stabilized samples, such as formalin fixed and paraffin-embedded (FFPE samples) and other samples that were treated with cross-linking fixatives such as glutaraldehyde.

Exemplary samples from which nucleic acids may be prepared include, but are not limited to, blood, swabs, body fluid, tissues including but not limited to, liver, spleen, kidney, lung, intestine, brain, heart, muscle, and pancreas, cell cultures, food samples, plant tissues or samples, as well as lysates, extracts, or materials and fractions obtained from the samples described above, or any cells, microorganisms and viruses that may be present on or in a sample, and the like.

Isolating target nucleic acids from a sample of interest may be performed by any method known in the art useful for nucleic acid isolation or purification. In addition, many kits for nucleic acid preparation are commercially available and may be used, including QIAamp DNA mini kit, QIAamp FFPE Tissue kit, and PAXgene DNA kit.

Alternatively, target nucleic acids useful as templates in single primer extension reactions disclosed herein may be obtained by modifying nucleic acids isolated from a sample. For example, in certain embodiments, target nucleic acids may be generated by (i) fragmenting DNA molecules (e.g., genomic DNA) to obtain double-stranded DNA fragments, and (ii) ligating a single-stranded common sequence to the 5' terminus of each stand of the double-stranded DNA fragments.

As indicated above, in certain embodiments, a single primer extension reaction may further comprise a boosting primer. A boosting primer may be included in a multiplex single primer extension reaction where most (i.e., at least about 90%, such as at least about 95%, about 98% or about 99%) or all of target nucleic acids each comprise a common sequence. The boosting primer comprises the common sequence or a portion thereof that is at least about 10 nucleotides in length, and is thus able to specifically anneal to a nucleic acid strand that comprises the complement of the common sequence and be extended using such a strand as a template in the single primer extension reaction. Thus, the presence of a boosting primer allows exponential amplification of target nucleic acids. However, because it comprises a sequence common to most or all of target nucleic acids, it does not affect specificity of the amplification. An exemplary use of a boosting primer is shown in FIG. 1 as described below in detail.

The initial concentration of a boosting primer may be from about 10 nM to about 10 mM, preferably from about 100 nM to about 2 mM, such as about 200 nM to about 600 nM.

Single primer extension reactions disclosed herein are performed in the presence of a MutS protein.

A MutS protein is a protein that binds to heteroduplex DNA containing mispaired or unpaired bases or a small DNA insertion/deletion loops of one to a few (e.g., 2, 3, 4, or 5) nucleotides and has sequence homology with *Escherichia coli* MutS protein (GenBank Accession No. GI: 127556, SEQ ID NO:1). It includes members of the MutS family described in Eisen, Nucleic Acids Research 26: 4291-4300, 1998, especially those in the MutS-I lineage.

In certain embodiments, the MutS protein is thermostable, that is, capable of at least partially maintaining (e.g., at least about 10%, about 20%, about 30%, about 40%, or about 50% of) its activity of binding to heteroduplex DNA containing mispaired or unpaired bases or a small DNA insertion/deletion loop of one to a few nucleotides at a temperature of at least about 50° C., such as at least about 60° C., about 70° C., about 75° C., about 80° C., about 85° C., about 90° C., or about 95° C. Such MutS protein may be of thermophilic or hyperthermophilic bacterial origin. Exemplary thermostable MutS proteins include MutS proteins from *Aquifex aeolicus*, *Aquifex pyrophilus*, *Thermotoga maritima*, *Thermus thermophilus* and *Thermus aquaticus* (see Fukui et al., Int. J. Mol. Sci. 14:6436-6453, 2013; U.S. Pat. No. 5,877,280; Biswas and Hsieh, Journal of Biological Chemistry 271:5040-5048, 1996; and Takamatsu et al., Nucleic Acid Research 24:640-647, 1996).

In certain preferred embodiments, the MutS protein is the MutS protein from *Thermus aquaticus* (herein referred to as "Taq MutS") (GenBank Accession No. GI: 1203807, SEQ ID NO:2). In related embodiments, the MutS protein is a MutS protein with a sequence identity of at least about 30%, such as at least about 40%, about 45%, about 50%, about 55%, about 60%, about 65%, about 70%, about 75%, about 80%, about 85%, about 90%, about 95%, about 96%, about 97%, about 98%, or about 99%, with Taq MutS.

In certain preferred embodiments, the MutS protein is the MutS protein from *Thermus thermophilus* (herein referred to as "Tth MutS") (GenBank Accession No. GI: 1871501, SEQ ID NO:3). In related embodiments, the MutS protein is a MutS protein having a sequence identity of at least about 30%, such as at least about 40%, about 45%, about 50%, about 55%, about 60%, about 65%, about 70%, about 75%, about 80%, about 85%, about 90%, about 95%, about 96%, about 97%, about 98%, or about 99%, with Tth MutS.

In certain preferred embodiments, the MutS protein is the MutS protein from *Aquifex aeolicus* (herein referred to as "Aae MutS") (GenBank Accession No. GI: 2983001, SEQ ID NO:4). In related embodiments, the MutS protein is a MutS protein having a sequence identity of at least about 30%, such as at least about 40%, about 45%, about 50%, about 55%, about 60%, about 65%, about 70%, about 75%, about 80%, about 85%, about 90%, about 95%, about 96%, about 97%, about 98%, or about 99%, with Aae MutS.

In certain preferred embodiments, the MutS protein is the MutS protein from *Aquifex pyrophilus* (herein referred to as "Apy MutS") (GenBank Accession No. GI: 1619907, SEQ ID NO:5). In related embodiments, the MutS protein is a MutS protein having a sequence identity of at least about 30%, such as at least about 40%, about 45%, about 50%, about 55%, about 60%, about 65%, about 70%, about 75%, about 80%, about 85%, about 90%, about 95%, about 96%, about 97%, about 98%, or about 99%, with Apy MutS.

In certain preferred embodiments, the MutS protein is the MutS protein from *Thermotoga maritima* (herein referred to as "Tma MutS") (GenBank Accession No. GI: 1619909, SEQ ID NO:6). In related embodiments, the MutS protein is a MutS protein having a sequence identity of at least about 30%, such as at least about 40%, about 45%, about 50%, about 55%, about 60%, about 65%, about 70%, about 75%, about 80%, about 85%, about 90%, about 95%, about 96%, about 97%, about 98%, or about 99%, with Tma MutS.

For determining sequence identity, the WU-BLAST-2 program (Altschul et al., Methods in Enzymology 266:460-480, 1996) is used. This program uses several search parameters, most of which are set to the default values. The adjustable parameters are set with the following values: overlap span=1, overlap fraction=0.125, word threshold (T)=11. The HSP S and HSP S2 parameters are dynamic values and are established by the program itself depending upon the composition of the particular sequence and composition of the particular database against which the sequence of interest is being searched. A percent sequence identity value is determined by the number of matching identical residues divided by the total number of residues of the "longer" sequence in the aligned region. The "longer" sequence is the one having the most actual residues in the aligned region (gaps introduced by WU-Blast-2 to maximize the alignment score are ignored).

MutS proteins may be obtained from an organism in which they occur in nature. Alternatively, they can be produced by recombinant methods.

A MutS protein may be present in a highly multiplex single primer extension reaction at a concentration of about 0.05 to about 10 uM, such as about 0.1 to about 1 uM, or about 0.2 to about 0.5 uM.

In addition to a MutS protein, single primer extension as disclosed herein is performed in the presence of Mg^{2+} at a minimum concentration of about 6 mM.

In certain embodiments, the concentration of Mg^{2+} is at least about 6, at least about 7, at least about 8, at least about 9, or at least about 10 mM and/or at most about 15, at most about 16, at most about 17, at most about 18, at most about 19, at most about 20 mM, including any combination of the above-noted minimum and maximum concentrations, such as about 6 mM to about 16 mM, about 6 mM to about 20 mM, about 8 mM to about 16 mM, about 8 mM to about 20 mM, about 10 mM to about 16 mM, or about 10 mM to about 20 mM.

Single primer extension reactions disclosed herein also comprise dNTPs. Typically, dNTPs are each present at an initial concentration (i.e., concentration at the start of the reaction) of about 40 uM to about 1000 uM in the reaction mixtures.

Single primer extension reactions disclosed herein also comprise a DNA polymerase for primer extensions. Preferably, DNA polymerases are thermostable. Exemplary DNA polymerases include Taq polymerase (from *Thermus aquaticus*), Tfi polymerase (from *Thermus filiformis*), Bst polymerase (from *Bacillus stearothermophilus*), Pfu polymerase (from *Pyrococcus furiosus*), Tth polymerase (from *Thermus thermophilus*), Pow polymerase (from *Pyrococcus woesei*), Tli polymerase (from *Thermococcus litoralis*), Ultima polymerase (from *Thermotoga maritima*), KOD polymerase (from *Thermococcus kodakaraensis*), Pol I and II polymerases (from *Pyrococcus abyssi*) and Pab (from *Pyrococcus abyssi*), Amplitaq Gold® DNA polymerase (Applied Biosciences), Stoffel fragment of Amplitaq® DNA Polymerase (Roche), KOD polymerase (EMD Biosciences), KOD Hot Start polymerase (EMD Biosciences), Deep

Vent™ DNA polymerase (New England Biolabs), Phusion polymerase (New England Biolabs), KlenTaq1 polymerase (DNA Polymerase Technology, Inc), KlenTaq Long Accuracy polymerase (DNA Polymerase Technology, Inc), Omni KlenTaq™ DNA polymerase (DNA Polymerase Technology, Inc), Omni KlenTaq™ LA DNA polymerase (DNA Polymerase Technology, Inc), Platinum® Taq DNA Polymerase (Invitrogen), Hemo KlenTaq™ (New England Biolabs), Platinum® Taq DNA Polymerase High Fidelity (Invitrogen), Platinum® Pfx (Invitrogen), Accuprime™ Pfx (Invitrogen), and Accuprime™ Taq DNA Polymerase High Fidelity (Invitrogen).

In embodiments where two or more primers are able to anneal to closely located regions of a target nucleic acid, the DNA polymerase preferably does not have strand displacement activity, flap endonuclease or 5'→3' exonuclease activity, with which the polymerase may destroy downstream primer extension products.

In certain preferred embodiments, the DNA polymerase does not have 3'→5" exonuclease activity.

As disclosed above, the combination of a relatively high concentration of Mg^{2+} and a MutS improves primer annealing and extension specificity in a highly multiplex single primer extension reaction. In certain embodiments, the primer annealing and extension specificity is at least about 80%, about 85%, about 90% or about 95%.

The specificity of primer annealing and extension in highly multiplex single primer extension can be measured, after next generation sequencing, by the percentage of reads mapped to the targeted loci.

Also as disclosed above, the combination of a relatively high concentration of Mg^{2+} and a MutS also improves primer annealing and extension efficiency in a highly multiplex primer extension reaction. In certain embodiments, the average primer annealing and extension efficiency of all of the primers other than a boosting primer (if present) in a highly multiplex single primer extension reaction is at least about 5%, about 10%, about 15%, about 20%, or about 25%.

The efficiency of primer annealing of an individual enrichment primer in a highly multiplex single primer extension reaction can be measured by measuring the number of DNA copies generated by the individual enrichment primer after primer extension, such as by a qPCR assay. For example, in certain embodiments where each enrichment primer comprises at its 5' terminus a universal sequence that is incapable of specifically annealing to a target nucleic acid, one qPCR primer can be designed to specifically anneal to the universal sequence and another qPCR primer can be designed to have a sequence identical to a portion of the target nucleic acid 5' to the region where the enrichment primer anneals. Such a qPCR assay measures the number of new DNA molecules generated by primer extension. The ratio of new DNA molecules to the input DNA amount represents the percentage of input DNA being annealed and extended by the enrichment primer.

In certain embodiments, the primer annealing and extension specificity is at least about 80%, about 85%, about 90% or about 95% and the average primer annealing and extension efficiency is at least about 5%, about 10%, about 15%, about 20%, or about 25%, including all possible combinations of the above-noted specificity and efficiency.

Products of single primer extension reactions disclosed herein may be further amplified and/or analyzed (e.g., sequenced). The further amplification may be performed using any appropriate methods known in the art.

Sequencing of extension products or amplified extension products may also be performed using any appropriate

methods known in the art. In certain embodiments, high throughput sequencing platforms known in the art may be used to sequence nucleic acids prepared from single primer extension products (see, Myllykangas et al., *Bioinformatics for High Throughput Sequencing*, Rodriguez-Ezpeleta et al. (eds.), Springer Science+Business Media, LLC, 2012, pages 11-25). Exemplary high throughput DNA sequencing systems include, but are not limited to, the GS FLX sequencing system originally developed by 454 Life Sciences and later acquired by Roche (Basel, Switzerland), Genome Analyzer developed by Solexa and later acquired by Illumina Inc. (San Diego, Calif.) (see, Bentley, *Curr Opin Genet Dev* 16:545-52, 2006; Bentley et al., *Nature* 456:53-59, 2008), the SOLiD sequence system by Life Technologies (Foster City, Calif.) (see, Smith et al., *Nucleic Acid Res* 38: e142, 2010; Valouev et al., *Genome Res* 18:1051-63, 2008), CGA developed by Complete Genomics and acquired by BGI (see, Drmanac et al., *Science* 327:78-81, 2010), PacBio RS sequencing technology developed by Pacific Biosciences (Menlo Park, Calif.) (see, Eid et al., *Science* 323: 133-8, 2009), and Ion Torrent developed by Life Technologies Corporation (see, U.S. Patent Application Publication Nos. 2009/0026082; 2010/0137143; and 2010/0282617).

The analysis of the extension products may be used in identifying structural variations (e.g., nucleotide mutations, deletions, insertions) in target nucleic acids when compared with corresponding reference sequences. In addition, sequencing extension products from cDNA is useful in detecting fusions and splicing variants. Furthermore, comparing sequences of extension products from bisulfite treated genomic DNA with those from untreated genomic DNA may be used in determining methylation patterns of the genomic DNA.

II. Using Single Primer Extension in NGS

In a related aspect, the present disclosure provides a method for using highly multiplex single primer extension reactions described above in high throughput nucleic acid sequencing. The method comprises:

a. providing a plurality of partially double-stranded target nucleic acids each of which comprises a double-stranded target nucleic acid sequence and a single-stranded common sequence at the 5' terminus of each strand of the target nucleic acid sequence, wherein the single-stranded common sequence is the same among the plurality of the target nucleic acids;

b. providing a plurality of primers capable of specifically annealing to portions of the target nucleic acid sequences;

c. extending the primers using the target nucleic acid sequences as templates in the presence of (1) a MutS protein and (2) Mg^{2+} at a minimum concentration of about 6 mM to obtain extension products; and

d. sequencing the extension products or their amplification products.

The partially double-stranded target nucleic acids may be prepared from any nucleic acids of interest (e.g., genomic DNA and bisulfite-treated genomic DNA) isolated from any samples, sources, or organisms as described above in Section I.

In certain embodiments, extending the primers in step c. is performed in the presence of a boosting primer comprising the common sequence or a portion thereof that is at least about 10 nucleotides in length.

In certain embodiments, step a. comprises:

a1. fragmenting double-stranded DNA molecules to generate double-stranded DNA fragments, and

a2. ligating the single-stranded common sequence to the 5' terminus of each strand of the double-stranded DNA fragment to provide the plurality of partially double-stranded nucleic acids.

Steps a1 and a2 may be performed using any appropriate methods known in the art. For example, step a2 may be performed using an adapter with the 5' portion in single-stranded formation and containing the common sequence and its 3' end portion forming a short double-stranded structure with a complementary oligonucleotide. The double-stranded 3' end can be ligated to the 5' terminus of each strand of the double-stranded DNA fragments in a regular double-strand DNA ligation reaction.

The single-stranded common sequence comprises at least about 10 nucleotides (e.g., at least about 11, 12, 13, 14, 15, 16, 17, 18, 19, or 20 nucleotides) so that a boosting primer may be designed to specifically anneal to the common sequence.

The plurality of primers in step b., single primer extensions in step c., and sequencing in step d. may be performed as described above in Section I.

In certain embodiments, the method also comprises: prior to step d., a further step of amplifying the extension products of step c. The amplification may be performed using a primer at least substantially complementary to the common sequence or a portion of the common sequence that is at least about 10 nucleotides in length and another primer capable of specifically annealing to a universal sequence at the 5' ends of the primers that have been extended in the single primer extension reaction.

An exemplary method is shown in FIG. 1. Briefly, in step 1, genomic DNA is first fragmented into smaller fragment (see step 1). In step 2, each fragment is then 5' ligated with a sequencing adapter ("NGSCounter") (an exemplary common sequence). In step 3, target-specific primers each having the same 5' universal sequence are annealed to the fragmented DNA and extended by DNA polymerase to the sequencing adapter ends. This process can be repeated multiple cycles to maximize the yield of single primer extension. An adapter primer ("boosting primer") can also be included to boost the yield of single primer extension. In step 4, the newly synthesized DNA strand, containing the 5' universal sequence and 3' adapter sequence, can be further amplified by a pair of an adapter primer (that comprises the boosting primer sequence or a portion thereof) and a universal primer (that comprises the universal sequence or a portion thereof) to prepare universal PCR library. In step 5, the library is quantified. In step 6, the nucleic acids in the library are sequenced.

III. Kits for Performing Single Primer Extensions

In another aspect, the present disclosure provides a kit for performing a highly multiplex single primer extension reaction, comprising: (1) a pool of at least 1,000 different primers, (2) a buffer concentrate comprising a MutS protein, dNTPs, and Mg^{2+} , wherein the concentration of Mg^{2+} in a highly multiplex single primer extension reaction mixture generated from diluting the buffer concentrate is at least about 6 mM, and (3) a DNA polymerase.

The pool of primers in the kit may contain about 1,000 to about 100,000 different primers. For example, a primer pool may contain at least about 1,000, 1,500, 2,000, 2,500, 3,000, 3,500, 4,000, 4,500, or 5,000 different primers, and/or at most about 10,000, 20,000, 30,000, 40,000, 50,000, 60,000, 70,000, 80,000, 90,000, or 100,000 different primers, including any combination of the above-noted minimum and maximum numbers of different primers such as about 1,000 to about 10,000, about 1,000 to about 50,000, about 1,000 to

about 100,000, about 5,000 to about 20,000, about 5,000 to about 50,000, or about 5,000 to about 100,000 different primers.

The primers in the pool may preferably be at the same concentration or a substantially similar concentration (i.e., the highest concentration of a primer is at most twice the lowest concentration of another primer).

Each primer in the pool may be about 0.01 nM to about 100 nM, preferably about 1 nM to about 40 nM, such as about 2 nM to about 20 nM.

In certain embodiments, each primer in the pool has a concentration of at least about 0.01, 0.02, 0.03, 0.04, 0.05, 0.06, 0.07, 0.08, 0.09, 0.1, 0.2, 0.3, 0.4, or 0.5 nM and/or at most about 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 30, 40, 50, 60, 70, 80, 90, or 100 nM, including any combination of the above-noted minimum and maximum concentrations.

The buffer concentrate in the kit comprises a MutS protein, dNTPs, and Mg^{2+} , wherein the concentration of Mg^{2+} in a highly multiplex single primer extension reaction mixture generated from diluting the buffer concentrate is at least about 6 mM, such as about 6 mM to about 20 mM.

In certain embodiments, the concentration of Mg^{2+} in a highly multiplex single primer extension reaction diluted from a buffer concentrate is at least about 6, at least about 7, at least about 8, at least about 9, or at least about 10 mM and/or at most about 15, at most about 16, at most about 17, at most about 18, at most about 19, at most about 20 mM, including any combination of the above-noted minimum and maximum concentrations, such as about 6 mM to about 16 mM, about 6 mM to about 20 mM, about 8 mM to about 16 mM, about 8 mM to about 20 mM, about 10 mM to about 16 mM, or about 10 mM to about 20 mM.

A MutS protein may be present in a highly multiplex single primer extension reaction diluted from a buffer concentrate at a concentration of about 0.05 to about 10 uM, such as about 0.1 to about 1 uM, or about 0.2 to about 0.5 uM.

dNTPs are typically each present at a concentration of about 40 to about 1000 uM in a highly multiplex single primer extension reaction diluted from a buffer concentrate.

Additional descriptions of various components of the kit, including primers, MutS proteins, dNTPs, DNA polymerases may be found above in section I.

The kit may further comprise one or more of the following additional components: a single-stranded sequence that may be ligated to the 5' terminus of each strand of double-stranded target nucleic acids as described above in Sections I and II, a boosting primer as described above in Sections I and II, a DNA ligase (e.g., T4 DNA ligase and *E. coli* DNA ligase), a ligation buffer, a primer pair for amplifying single primer extension products (e.g., those described above in Section II), and a sequencing primer.

In a related aspect, the present disclosure provides use of a kit disclosed herein in performing a highly multiplex primer extension reaction.

In another related aspect, the present disclosure provides use of a kit disclosed herein in nucleic acid sequencing.

EXAMPLES

Example 1

Effects of Mg^{2+} Concentrations on Primer Extension Efficiency and Specificity

Extension reactions of primers (each primer at a 2 nM concentration) were performed at Mg^{2+} concentrations from

2 mM to 20 mM. The primer extension reaction was set up in 1× miScript buffer supplemented with various amount of Mg^{2+} , 2.4 U HotStatTaq, 20 nM single target specific primer with 5' universal sequence and 20 ng genomic DNA. The reaction proceeded for 20 minutes at 95° C., 15 minutes at 55° C., 15 minutes at 65° C. and 7 minutes at 72° C. Specific qPCR was used to quantify the amount of extension products.

The results (FIG. 2) show that increasing the concentration of Mg^{2+} improved the primer annealing and extension efficiency.

DNA enriched using 12 mM Mg^{2+} and over 1300 primers at 2 nM each and DNA enriched using 4 mM Mg^{2+} and the same primers at 20 nM each were sequenced to determine primer annealing and extension specificity.

The results (FIG. 3) show that the higher Mg^{2+} concentration significantly reduced the specificity of primer annealing. While not wishing to be bound by any particular theory, the present inventors believe that under higher Mg^{2+} concentrations, non-specific binding of primers to DNA templates was stabilized along with increased specific annealing.

Example 2

Effects of MutS on Specificity of Primer Extensions Performed at Relatively High Mg^{2+} Concentrations

This example demonstrates that MutS improved specificity of highly multiplex targeted enrichment performed at a relatively high Mg^{2+} concentration.

Methods

Fragmenting gDNA (Sonication Method)

1 ug genomic DNA (gDNA) was diluted in 100 ul elution buffer and sonicated on BIORUPTOR300 (power level low, 10 cycles of 30 seconds on and 30 seconds off). The resulting DNA fragments were purified using QIAGEN's MinElute Column.

Library Construction

A sequencing library was prepared using the above-generated gDNA fragments according to the End Repair, A-Addition, Adapter Ligation and Clean up and Size Selection of DNA protocols as described in QIAGEN GeneRead Library Prep for Illumina Kit. The library was quantified and diluted to 20 ng/ul.

Primer Extension

The primer extension reactions were set up as follows:

Primer Extension	20 nM Primer Setup Volume (ul)	2 nM Primer setup Volume (ul)
GeneRead V2 buffer (5x)	8	8
dNTP (2 mM each)	4	4
DNA library (from previous step, 20 ng/ul)	1	1
Taq	1.6	1.6
SPE Primer pool (76 nM)	10.5 20 nM final	1.05 2 nM final
Illumina_F primer (10 uM)	1.6	1.6
MgCl ₂ (80 mM)	0	4
Taq MutS	0	3
H ₂ O	13.3	15.75
Total	40	40

The final Mg^{2+} concentrations of 20 nM primer setup and 2 nM primer setup were 4 and 12 mM, respectively.

15

Primer extension reactions were performed in a thermal cycler as follows: 95° C. for 15 minutes, 5 cycles of 95° C. for 15 seconds and 68° C. for 10 minutes, 72° C. for 5 minutes, and 4° C. hold.

60 ul H₂O was added to each of the primer extension mixtures to bring the volume to 100 ul. The extension products were purified by 2 rounds Ampure beads purification (1.2× beads ratio) and eluted in 15 ul Buffer EB. Universal PCR

Primer extension products were further amplified in a universal PCR reaction as follows:

Universal qPCR	Volume (ul)
GeneRead V2 buffer (5x)	5
dNTP (2 mM each)	2.5
Template (from previous step)	11.5
Taq	1
Universal primer (IL2) (4 uM)	2.5
IL1_ID(#)-RS2 primer (4 uM)	2.5
Total	25

The reaction was performed in a thermal cycler as follows: 95° C. for 15 minutes, 26 cycles of 95° C. for 15 seconds and 65° C. for 2 minutes, and 4° C. hold.

The amplification products were purified via one round Ampure beads purification (1.2× beads ratio) and eluted in 20 ul buffer EB. 5 ul of the purified amplification products were checked with a 1.2% Agarose gel.

Library Quantification and Sequencing

Quantification of the amplified library and sequencing were performed according to the manufacturer's instructions.

Results

The results (FIG. 4) show that adding Taq MutS in primer extension reactions together with 12 mM Mg²⁺ significantly improve primer annealing and extension specificity at no

16

cost of primer annealing and extension efficiency as indicated by consistently good sequence uniformity.

T50 is an indicator of sequence uniformity. It is defined as the percentage of total sequence throughput used by regions sequenced at below 50% of the mean read depth of all target regions. In a perfect uniform scenario, the T50 value is 50.

While not wishing to be bound by any particular theory, the present inventors hypothesize that non-specific primer-template hybrids formed during primer extension performed at a high Mg²⁺ concentration presumably contain mismatches and bulges in the middle of the primer sequences, but perfectly anneal at the 3' termini of the primers so that DNA polymerase can extend from such termini. MutS specifically binds to those bulges created by sequence mismatches and thus block a DNA polymerase from binding to and extending the non-specific primer-template duplex. The specific primer-template duplex would have complementary sequences and no bulge in the structure so that MutS would not bind tightly. Thus, extension of primers that specifically bind to their templates will not be affected.

The various embodiments described above can be combined to provide further embodiments. All of the U.S. patents, U.S. patent application publications, U.S. patent applications, foreign patents, foreign patent applications and non-patent publications referred to in this specification and/or listed in the Application Data Sheet are incorporated herein by reference, in their entirety. Aspects of the embodiments can be modified, if necessary to employ concepts of the various patents, applications and publications to provide yet further embodiments.

These and other changes can be made to the embodiments in light of the above-detailed description. In general, in the following claims, the terms used should not be construed to limit the claims to the specific embodiments disclosed in the specification and the claims, but should be construed to include all possible embodiments along with the full scope of equivalents to which such claims are entitled. Accordingly, the claims are not limited by the disclosure.

SEQUENCE LISTING

<160> NUMBER OF SEQ ID NOS: 6

<210> SEQ ID NO 1

<211> LENGTH: 853

<212> TYPE: PRT

<213> ORGANISM: Escherichia coli

<400> SEQUENCE: 1

Met Ser Ala Ile Glu Asn Phe Asp Ala His Thr Pro Met Met Gln Gln
1 5 10 15

Tyr Leu Arg Leu Lys Ala Gln His Pro Glu Ile Leu Leu Phe Tyr Arg
20 25 30

Met Gly Asp Phe Tyr Glu Leu Phe Tyr Asp Asp Ala Lys Arg Ala Ser
35 40 45

Gln Leu Leu Asp Ile Ser Leu Thr Lys Arg Gly Ala Ser Ala Gly Glu
50 55 60

Pro Ile Pro Met Ala Gly Ile Pro Tyr His Ala Val Glu Asn Tyr Leu
65 70 75 80

Ala Lys Leu Val Asn Gln Gly Glu Ser Val Ala Ile Cys Glu Gln Ile
85 90 95

Gly Asp Pro Ala Thr Ser Lys Gly Pro Val Glu Arg Lys Val Val Arg
100 105 110

-continued

Ile Val Thr Pro Gly Thr Ile Ser Asp Glu Ala Leu Leu Gln Glu Arg	
115	120 125
Gln Asp Asn Leu Leu Ala Ala Ile Trp Gln Asp Ser Lys Gly Phe Gly	
130	135 140
Tyr Ala Thr Leu Asp Ile Ser Ser Gly Arg Phe Arg Leu Ser Glu Pro	
145	150 155 160
Ala Asp Arg Glu Thr Met Ala Ala Glu Leu Gln Arg Thr Asn Pro Ala	
	165 170 175
Glu Leu Leu Tyr Ala Glu Asp Phe Ala Glu Met Ser Leu Ile Glu Gly	
	180 185 190
Arg Arg Gly Leu Arg Arg Arg Pro Leu Trp Glu Phe Glu Ile Asp Thr	
	195 200 205
Ala Arg Gln Gln Leu Asn Leu Gln Phe Gly Thr Arg Asp Leu Val Gly	
	210 215 220
Phe Gly Val Glu Asn Ala Pro Arg Gly Leu Cys Ala Ala Gly Cys Leu	
225	230 235 240
Leu Gln Tyr Ala Lys Asp Thr Gln Arg Thr Thr Leu Pro His Ile Arg	
	245 250 255
Ser Ile Thr Met Glu Arg Glu Gln Asp Ser Ile Ile Met Asp Ala Ala	
	260 265 270
Thr Arg Arg Asn Leu Glu Ile Thr Gln Asn Leu Ala Gly Gly Ala Glu	
	275 280 285
Asn Thr Leu Ala Ser Val Leu Asp Cys Thr Val Thr Pro Met Gly Ser	
	290 295 300
Arg Met Leu Lys Arg Trp Leu His Met Pro Val Arg Asp Thr Arg Val	
305	310 315 320
Leu Leu Glu Arg Gln Gln Thr Ile Gly Ala Leu Gln Asp Phe Thr Ala	
	325 330 335
Gly Leu Gln Pro Val Leu Arg Gln Val Gly Asp Leu Glu Arg Ile Leu	
	340 345 350
Ala Arg Leu Ala Leu Arg Thr Ala Arg Pro Arg Asp Leu Ala Arg Met	
	355 360 365
Arg His Ala Phe Gln Gln Leu Pro Glu Leu Arg Ala Gln Leu Glu Thr	
	370 375 380
Val Asp Ser Ala Pro Val Gln Ala Leu Arg Glu Lys Met Gly Glu Phe	
385	390 395 400
Ala Glu Leu Arg Asp Leu Leu Glu Arg Ala Ile Ile Asp Thr Pro Pro	
	405 410 415
Val Leu Val Arg Asp Gly Gly Val Ile Ala Ser Gly Tyr Asn Glu Glu	
	420 425 430
Leu Asp Glu Trp Arg Ala Leu Ala Asp Gly Ala Thr Asp Tyr Leu Glu	
	435 440 445
Arg Leu Glu Val Arg Glu Arg Glu Arg Thr Gly Leu Asp Thr Leu Lys	
	450 455 460
Val Gly Phe Asn Ala Val His Gly Tyr Tyr Ile Gln Ile Ser Arg Gly	
465	470 475 480
Gln Ser His Leu Ala Pro Ile Asn Tyr Met Arg Arg Gln Thr Leu Lys	
	485 490 495
Asn Ala Glu Arg Tyr Ile Ile Pro Glu Leu Lys Glu Tyr Glu Asp Lys	
	500 505 510
Val Leu Thr Ser Lys Gly Lys Ala Leu Ala Leu Glu Lys Gln Leu Tyr	
	515 520 525

-continued

Glu Glu Leu Phe Asp Leu Leu Leu Pro His Leu Glu Ala Leu Gln Gln
 530 535 540
 Ser Ala Ser Ala Leu Ala Glu Leu Asp Val Leu Val Asn Leu Ala Glu
 545 550 555 560
 Arg Ala Tyr Thr Leu Asn Tyr Thr Cys Pro Thr Phe Ile Asp Lys Pro
 565 570 575
 Gly Ile Arg Ile Thr Glu Gly Arg His Pro Val Val Glu Gln Val Leu
 580 585 590
 Asn Glu Pro Phe Ile Ala Asn Pro Leu Asn Leu Ser Pro Gln Arg Arg
 595 600 605
 Met Leu Ile Ile Thr Gly Pro Asn Met Gly Gly Lys Ser Thr Tyr Met
 610 615 620
 Arg Gln Thr Ala Leu Ile Ala Leu Met Ala Tyr Ile Gly Ser Tyr Val
 625 630 635 640
 Pro Ala Gln Lys Val Glu Ile Gly Pro Ile Asp Arg Ile Phe Thr Arg
 645 650 655
 Val Gly Ala Ala Asp Asp Leu Ala Ser Gly Arg Ser Thr Phe Met Val
 660 665 670
 Glu Met Thr Glu Thr Ala Asn Ile Leu His Asn Ala Thr Glu Tyr Ser
 675 680 685
 Leu Val Leu Met Asp Glu Ile Gly Arg Gly Thr Ser Thr Tyr Asp Gly
 690 695 700
 Leu Ser Leu Ala Trp Ala Cys Ala Glu Asn Leu Ala Asn Lys Ile Lys
 705 710 715 720
 Ala Leu Thr Leu Phe Ala Thr His Tyr Phe Glu Leu Thr Gln Leu Pro
 725 730 735
 Glu Lys Met Glu Gly Val Ala Asn Val His Leu Asp Ala Leu Glu His
 740 745 750
 Gly Asp Thr Ile Ala Phe Met His Ser Val Gln Asp Gly Ala Ala Ser
 755 760 765
 Lys Ser Tyr Gly Leu Ala Val Ala Ala Leu Ala Gly Val Pro Lys Glu
 770 775 780
 Val Ile Lys Arg Ala Arg Gln Lys Leu Arg Glu Leu Glu Ser Ile Ser
 785 790 795 800
 Pro Asn Ala Ala Ala Thr Gln Val Asp Gly Thr Gln Met Ser Leu Leu
 805 810 815
 Ser Val Pro Glu Glu Thr Ser Pro Ala Val Glu Ala Leu Glu Asn Leu
 820 825 830
 Asp Pro Asp Ser Leu Thr Pro Arg Gln Ala Leu Glu Trp Ile Tyr Arg
 835 840 845
 Leu Lys Ser Leu Val
 850

<210> SEQ ID NO 2

<211> LENGTH: 811

<212> TYPE: PRT

<213> ORGANISM: Thermus aquaticus

<400> SEQUENCE: 2

Met Glu Gly Met Leu Lys Gly Glu Gly Pro Gly Pro Leu Pro Pro Leu
 1 5 10 15
 Leu Gln Gln Tyr Val Glu Leu Arg Asp Gln Tyr Pro Asp Tyr Leu Leu
 20 25 30
 Leu Phe Gln Val Gly Asp Phe Tyr Glu Cys Phe Gly Glu Asp Ala Glu
 35 40 45

-continued

Arg Leu Ala Arg Ala Leu Gly Leu Val Leu Thr His Lys Thr Ser Lys
 50 55 60
 Asp Phe Thr Thr Pro Met Ala Gly Ile Pro Leu Arg Ala Phe Glu Ala
 65 70 75 80
 Tyr Ala Glu Arg Leu Leu Lys Met Gly Phe Arg Leu Ala Val Ala Asp
 85 90 95
 Gln Val Glu Pro Ala Glu Glu Ala Glu Gly Leu Val Arg Arg Glu Val
 100 105 110
 Thr Gln Leu Leu Thr Pro Gly Thr Leu Leu Gln Glu Ser Leu Leu Pro
 115 120 125
 Arg Glu Ala Asn Tyr Leu Ala Ala Ile Ala Thr Gly Asp Gly Trp Gly
 130 135 140
 Leu Ala Phe Leu Asp Val Ser Thr Gly Glu Phe Lys Gly Thr Val Leu
 145 150 155 160
 Lys Ser Lys Ser Ala Leu Tyr Asp Glu Leu Phe Arg His Arg Pro Ala
 165 170 175
 Glu Val Leu Leu Ala Pro Glu Leu Leu Glu Asn Gly Ala Phe Leu Asp
 180 185 190
 Glu Phe Arg Lys Arg Phe Pro Val Met Leu Ser Glu Ala Pro Phe Glu
 195 200 205
 Pro Glu Gly Glu Gly Pro Leu Ala Leu Arg Arg Ala Arg Gly Ala Leu
 210 215 220
 Leu Ala Tyr Ala Gln Arg Thr Gln Gly Gly Ala Leu Ser Leu Gln Pro
 225 230 235 240
 Phe Arg Phe Tyr Asp Pro Gly Ala Phe Met Arg Leu Pro Glu Ala Thr
 245 250 255
 Leu Arg Ala Leu Glu Val Phe Glu Pro Leu Arg Gly Gln Asp Thr Leu
 260 265 270
 Phe Ser Val Leu Asp Glu Thr Arg Thr Ala Pro Gly Arg Arg Leu Leu
 275 280 285
 Gln Ser Trp Leu Arg His Pro Leu Leu Asp Arg Gly Pro Leu Glu Ala
 290 295 300
 Arg Leu Asp Arg Val Glu Gly Phe Val Arg Glu Gly Ala Leu Arg Glu
 305 310 315 320
 Gly Val Arg Arg Leu Leu Tyr Arg Leu Ala Asp Leu Glu Arg Leu Ala
 325 330 335
 Thr Arg Leu Glu Leu Gly Arg Ala Ser Pro Lys Asp Leu Gly Ala Leu
 340 345 350
 Arg Arg Ser Leu Gln Ile Leu Pro Glu Leu Arg Ala Leu Leu Gly Glu
 355 360 365
 Glu Val Gly Leu Pro Asp Leu Ser Pro Leu Lys Glu Glu Leu Glu Ala
 370 375 380
 Ala Leu Val Glu Asp Pro Pro Leu Lys Val Ser Glu Gly Gly Leu Ile
 385 390 395 400
 Arg Glu Gly Tyr Asp Pro Asp Leu Asp Ala Leu Arg Ala Ala His Arg
 405 410 415
 Glu Gly Val Ala Tyr Phe Leu Glu Leu Glu Glu Arg Glu Arg Glu Arg
 420 425 430
 Thr Gly Ile Pro Thr Leu Lys Val Gly Tyr Asn Ala Val Phe Gly Tyr
 435 440 445
 Tyr Leu Glu Val Thr Arg Pro Tyr Tyr Glu Arg Val Pro Lys Glu Tyr
 450 455 460

-continued

Arg Pro Val Gln Thr Leu Lys Asp Arg Gln Arg Tyr Thr Leu Pro Glu
 465 470 475 480
 Met Lys Glu Lys Glu Arg Glu Val Tyr Arg Leu Glu Ala Leu Ile Arg
 485 490 495
 Arg Arg Glu Glu Glu Val Phe Leu Glu Val Arg Glu Arg Ala Lys Arg
 500 505 510
 Gln Ala Glu Ala Leu Arg Glu Ala Ala Arg Ile Leu Ala Glu Leu Asp
 515 520 525
 Val Tyr Ala Ala Leu Ala Glu Val Ala Val Arg Tyr Gly Tyr Val Arg
 530 535 540
 Pro Arg Phe Gly Asp Arg Leu Gln Ile Arg Ala Gly Arg His Pro Val
 545 550 555 560
 Val Glu Arg Arg Thr Glu Phe Val Pro Asn Asp Leu Glu Met Ala His
 565 570 575
 Glu Leu Val Leu Ile Thr Gly Pro Asn Met Ala Gly Lys Ser Thr Phe
 580 585 590
 Leu Arg Gln Thr Ala Leu Ile Ala Leu Leu Ala Gln Val Gly Ser Phe
 595 600 605
 Val Pro Ala Glu Glu Ala His Leu Pro Leu Phe Asp Gly Ile Tyr Thr
 610 615 620
 Arg Ile Gly Ala Ser Asp Asp Leu Ala Gly Gly Lys Ser Thr Phe Met
 625 630 635 640
 Val Glu Met Glu Glu Val Ala Leu Ile Leu Lys Glu Ala Thr Glu Asn
 645 650 655
 Ser Leu Val Leu Leu Asp Glu Val Gly Arg Gly Thr Ser Ser Leu Asp
 660 665 670
 Gly Val Ala Ile Ala Thr Ala Val Ala Glu Ala Leu His Glu Arg Arg
 675 680 685
 Ala Tyr Thr Leu Phe Ala Thr His Tyr Phe Glu Leu Thr Ala Leu Gly
 690 695 700
 Leu Pro Arg Leu Lys Asn Leu His Val Ala Ala Arg Glu Glu Ala Gly
 705 710 715 720
 Gly Leu Val Phe Tyr His Gln Val Leu Pro Gly Pro Ala Ser Lys Ser
 725 730 735
 Tyr Gly Val Glu Val Ala Ala Met Ala Gly Leu Pro Lys Glu Val Val
 740 745 750
 Ala Arg Ala Arg Ala Leu Leu Gln Ala Met Ala Ala Arg Arg Glu Gly
 755 760 765
 Ala Leu Asp Ala Val Leu Glu Arg Leu Leu Ala Leu Asp Pro Asp Arg
 770 775 780
 Leu Thr Pro Leu Glu Ala Leu Arg Leu Leu Gln Glu Leu Lys Ala Leu
 785 790 795 800
 Ala Leu Gly Ala Pro Leu Asp Thr Met Lys Gly
 805 810

<210> SEQ ID NO 3

<211> LENGTH: 819

<212> TYPE: PRT

<213> ORGANISM: Thermus thermophilus

<400> SEQUENCE: 3

Met Gly Gly Tyr Gly Gly Val Lys Met Glu Gly Met Leu Lys Gly Glu
 1 5 10 15
 Gly Pro Gly Pro Leu Pro Pro Leu Leu Gln Gln Tyr Val Glu Leu Arg
 20 25 30

-continued

Asp Arg Tyr Pro Asp Tyr Leu Leu Leu Phe Gln Val Gly Asp Phe Tyr
 35 40 45
 Glu Cys Phe Gly Glu Asp Ala Glu Arg Leu Ala Arg Ala Leu Gly Leu
 50 55 60
 Val Leu Thr His Lys Thr Ser Lys Asp Phe Thr Thr Pro Met Ala Gly
 65 70 75 80
 Ile Pro Ile Arg Ala Phe Asp Ala Tyr Ala Glu Arg Leu Leu Lys Met
 85 90 95
 Gly Phe Arg Leu Ala Val Ala Asp Gln Val Glu Pro Ala Glu Glu Ala
 100 105 110
 Glu Gly Leu Val Arg Arg Glu Val Thr Gln Leu Leu Thr Pro Gly Thr
 115 120 125
 Leu Thr Gln Glu Ala Leu Leu Pro Arg Glu Ala Asn Tyr Leu Ala Ala
 130 135 140
 Ile Ala Thr Gly Asp Gly Trp Gly Leu Ala Phe Leu Asp Val Ser Thr
 145 150 155 160
 Gly Glu Phe Lys Gly Thr Leu Leu Lys Ser Lys Ser Ala Leu Tyr Asp
 165 170 175
 Glu Leu Phe Arg His Arg Pro Ala Glu Val Leu Leu Ala Pro Glu Leu
 180 185 190
 Arg Glu Asn Glu Ala Phe Val Ala Glu Phe Arg Lys Arg Phe Pro Val
 195 200 205
 Met Leu Ser Glu Ala Pro Phe Glu Pro Gln Gly Glu Gly Pro Leu Ala
 210 215 220
 Leu Arg Arg Ala Gln Gly Ala Leu Leu Ala Tyr Ala Arg Ala Thr Gln
 225 230 235 240
 Gly Gly Ala Leu Ser Val Arg Pro Phe Arg Leu Tyr Asp Pro Gly Ala
 245 250 255
 Phe Val Arg Leu Pro Glu Ala Ser Leu Lys Ala Leu Glu Val Phe Glu
 260 265 270
 Pro Leu Arg Gly Gln Asp Thr Leu Phe Gly Val Leu Asp Glu Thr Arg
 275 280 285
 Thr Ala Pro Gly Arg Arg Leu Leu Gln Ala Trp Leu Arg His Pro Leu
 290 295 300
 Leu Glu Arg Gly Pro Leu Glu Ala Arg Leu Asp Arg Val Glu Arg Phe
 305 310 315 320
 Val Arg Glu Gly Ala Leu Arg Glu Gly Val Arg Arg Leu Leu Phe Arg
 325 330 335
 Leu Ala Asp Leu Glu Arg Leu Ala Thr Arg Leu Glu Leu Ser Arg Ala
 340 345 350
 Ser Pro Arg Asp Leu Ala Ala Leu Arg Arg Ser Leu Glu Ile Leu Pro
 355 360 365
 Glu Leu Lys Gly Leu Leu Gly Glu Glu Val Gly Leu Pro Asp Leu Ser
 370 375 380
 Gly Leu Leu Glu Glu Leu Arg Ala Ala Leu Val Glu Asp Pro Pro Leu
 385 390 395 400
 Lys Val Ser Glu Gly Gly Leu Ile Arg Glu Gly Tyr Asp Pro Asp Leu
 405 410 415
 Asp Ala Leu Arg Arg Ala His Ala Glu Gly Val Ala Tyr Phe Leu Asp
 420 425 430
 Leu Glu Ala Arg Glu Lys Glu Arg Thr Gly Ile Pro Thr Leu Lys Val
 435 440 445

-continued

Gly	Tyr	Asn	Ala	Val	Phe	Gly	Tyr	Tyr	Leu	Glu	Val	Thr	Arg	Pro	Tyr
450						455					460				
Tyr	Glu	Lys	Val	Pro	Gln	Glu	Tyr	Arg	Pro	Val	Gln	Thr	Leu	Lys	Asp
465					470					475					480
Arg	Gln	Arg	Tyr	Thr	Leu	Pro	Glu	Met	Lys	Glu	Arg	Glu	Arg	Glu	Leu
				485					490						495
Tyr	Arg	Leu	Glu	Ala	Leu	Ile	Lys	Arg	Arg	Glu	Glu	Glu	Val	Phe	Leu
		500						505					510		
Ala	Leu	Arg	Glu	Arg	Ala	Arg	Lys	Glu	Ala	Glu	Ala	Leu	Arg	Glu	Ala
		515					520						525		
Ala	Arg	Ile	Leu	Ala	Glu	Leu	Asp	Val	Tyr	Ala	Ala	Leu	Ala	Glu	Val
		530				535					540				
Ala	Val	Arg	His	Gly	Tyr	Thr	Arg	Pro	Arg	Phe	Gly	Glu	Arg	Leu	Arg
545					550					555					560
Ile	Arg	Ala	Gly	Arg	His	Pro	Val	Val	Glu	Arg	Arg	Thr	Ala	Phe	Val
			565						570					575	
Pro	Asn	Asp	Leu	Glu	Met	Ala	His	Glu	Leu	Val	Leu	Val	Thr	Gly	Pro
			580					585					590		
Asn	Met	Ala	Gly	Lys	Ser	Thr	Phe	Leu	Arg	Gln	Thr	Ala	Leu	Ile	Ala
		595					600					605			
Leu	Leu	Ala	Gln	Ile	Gly	Ser	Phe	Val	Pro	Ala	Glu	Glu	Ala	Glu	Leu
		610				615					620				
Pro	Leu	Phe	Asp	Gly	Ile	Tyr	Thr	Arg	Ile	Gly	Ala	Ser	Asp	Asp	Leu
625					630					635					640
Ala	Gly	Gly	Lys	Ser	Thr	Phe	Met	Val	Glu	Met	Glu	Glu	Val	Ala	Leu
			645						650					655	
Val	Leu	Lys	Glu	Ala	Thr	Glu	Arg	Ser	Leu	Val	Leu	Leu	Asp	Glu	Val
			660					665					670		
Gly	Arg	Gly	Thr	Ser	Ser	Leu	Asp	Gly	Val	Ala	Ile	Ala	Thr	Ala	Leu
		675					680					685			
Ala	Glu	Ala	Leu	His	Glu	Arg	Arg	Cys	Tyr	Thr	Leu	Phe	Ala	Thr	His
		690				695					700				
Tyr	Phe	Glu	Leu	Thr	Ala	Leu	Ala	Leu	Pro	Arg	Leu	Lys	Asn	Leu	His
705					710					715					720
Val	Ala	Ala	Lys	Glu	Glu	Glu	Gly	Gly	Leu	Val	Phe	Tyr	His	Gln	Val
			725						730					735	
Leu	Pro	Gly	Pro	Ala	Ser	Lys	Ser	Tyr	Gly	Val	Glu	Val	Ala	Glu	Met
			740					745					750		
Ala	Gly	Leu	Pro	Lys	Glu	Val	Val	Glu	Arg	Ala	Arg	Ala	Leu	Leu	Ser
		755					760					765			
Ala	Met	Ala	Ala	Arg	Arg	Glu	Gly	Ala	Leu	Glu	Glu	Val	Leu	Glu	Arg
		770				775					780				
Leu	Leu	Ala	Leu	Asp	Pro	Asp	Arg	Leu	Thr	Pro	Leu	Glu	Ala	Leu	Arg
785					790					795					800
Phe	Leu	His	Glu	Leu	Lys	Ala	Leu	Ala	Leu	Gly	Leu	Pro	Leu	Gly	Ser
			805					810						815	

Met Lys Gly

<210> SEQ ID NO 4
 <211> LENGTH: 859
 <212> TYPE: PRT
 <213> ORGANISM: Aquifex aeolicus

<400> SEQUENCE: 4

-continued

Met	Glu	Lys	Ser	Glu	Lys	Glu	Leu	Thr	Pro	Met	Leu	Ser	Gln	Tyr	His	
1				5					10					15		
Tyr	Phe	Lys	Asn	Gln	Tyr	Pro	Asp	Cys	Leu	Leu	Leu	Phe	Arg	Leu	Gly	
			20				25						30			
Asp	Phe	Tyr	Glu	Leu	Phe	Tyr	Glu	Asp	Ala	Tyr	Ile	Gly	Ser	Lys	Glu	
			35				40					45				
Leu	Gly	Leu	Val	Leu	Thr	Ser	Arg	Pro	Ala	Gly	Lys	Gly	Lys	Glu	Arg	
	50					55					60					
Ile	Pro	Met	Cys	Gly	Val	Pro	Tyr	His	Ser	Ala	Asn	Ser	Tyr	Ile	Ala	
65					70					75					80	
Lys	Leu	Val	Asn	Lys	Gly	Tyr	Lys	Val	Ala	Ile	Cys	Glu	Gln	Val	Glu	
			85						90					95		
Asp	Pro	Ser	Lys	Ala	Lys	Gly	Ile	Val	Lys	Arg	Glu	Val	Val	Arg	Val	
			100					105					110			
Ile	Thr	Pro	Gly	Thr	Phe	Phe	Glu	Arg	Asp	Thr	Gly	Gly	Leu	Ala	Ser	
	115						120					125				
Leu	Tyr	Lys	Lys	Gly	Asn	His	Tyr	Tyr	Val	Gly	Tyr	Leu	Asn	Leu	Ala	
	130					135					140					
Val	Gly	Glu	Phe	Leu	Gly	Ala	Lys	Val	Lys	Ile	Glu	Glu	Leu	Leu	Asp	
145					150					155					160	
Leu	Leu	Ser	Lys	Leu	Asn	Ile	Lys	Glu	Ile	Leu	Val	Lys	Lys	Gly	Glu	
			165					170						175		
Lys	Leu	Pro	Glu	Glu	Leu	Glu	Lys	Val	Leu	Lys	Val	Tyr	Val	Ser	Glu	
			180					185					190			
Leu	Glu	Glu	Glu	Phe	Phe	Glu	Glu	Gly	Ser	Glu	Glu	Ile	Leu	Lys	Asp	
	195						200					205				
Phe	Gly	Val	Leu	Ser	Leu	Gln	Ala	Phe	Gly	Phe	Glu	Glu	Asp	Thr	Tyr	
	210					215					220					
Ser	Leu	Pro	Leu	Gly	Ala	Val	Tyr	Lys	Tyr	Ala	Lys	Thr	Thr	Gln	Lys	
225				230						235					240	
Gly	Tyr	Thr	Pro	Leu	Ile	Pro	Arg	Pro	Lys	Pro	Tyr	Arg	Asp	Glu	Gly	
			245						250					255		
Phe	Val	Arg	Leu	Asp	Ile	Lys	Ala	Ile	Lys	Gly	Leu	Glu	Ile	Leu	Glu	
	260						265						270			
Ser	Leu	Glu	Gly	Arg	Lys	Asp	Ile	Ser	Leu	Phe	Lys	Val	Ile	Asp	Arg	
	275						280					285				
Thr	Leu	Thr	Gly	Met	Gly	Arg	Arg	Arg	Leu	Lys	Phe	Arg	Leu	Leu	Ser	
	290					295					300					
Pro	Phe	Arg	Ser	Arg	Glu	Lys	Ile	Glu	Arg	Ile	Gln	Glu	Gly	Val	Gln	
305					310					315					320	
Glu	Leu	Lys	Glu	Asn	Arg	Glu	Ala	Leu	Leu	Lys	Ile	Arg	Gln	Ile	Leu	
			325						330					335		
Glu	Gly	Met	Ala	Asp	Leu	Glu	Arg	Leu	Val	Ser	Lys	Ile	Ser	Ser	Asn	
			340					345				350				
Met	Ala	Thr	Pro	Arg	Glu	Leu	Val	Tyr	Leu	Lys	Asn	Ser	Leu	Lys	Lys	
	355						360					365				
Val	Glu	Glu	Leu	Arg	Leu	Leu	Leu	Leu	Glu	Leu	Lys	Ala	Pro	Ile	Phe	
	370					375					380					
Lys	Glu	Ile	Leu	Gln	Asn	Phe	Glu	Asp	Thr	Lys	Lys	Ile	Ile	Asn	Asp	
385					390					395					400	
Ile	Glu	Lys	Thr	Leu	Val	Glu	Asp	Pro	Pro	Leu	His	Val	Lys	Glu	Gly	
			405						410					415		
Gly	Leu	Ile	Arg	Glu	Gly	Val	Asn	Ala	Tyr	Leu	Asp	Glu	Leu	Arg	Phe	

-continued

420						425						430					
Ile	Arg	Asp	Asn	Ala	Glu	Thr	Tyr	Leu	Arg	Glu	Tyr	Glu	Lys	Lys	Leu		
		435					440					445					
Arg	Gln	Glu	Thr	Gly	Ile	Gln	Ser	Leu	Lys	Ile	Gly	Tyr	Asn	Lys	Val		
	450					455					460						
Met	Gly	Tyr	Tyr	Ile	Glu	Val	Thr	Lys	Pro	Asn	Leu	Lys	Tyr	Val	Pro		
465					470					475					480		
Ser	Tyr	Phe	Arg	Arg	Arg	Gln	Thr	Leu	Ser	Asn	Ser	Glu	Arg	Phe	Thr		
				485					490					495			
Thr	Glu	Glu	Leu	Gln	Arg	Leu	Glu	Glu	Lys	Ile	Leu	Ser	Ala	Gln	Thr		
			500					505					510				
Arg	Ile	Asn	Asp	Leu	Glu	Tyr	Glu	Leu	Tyr	Lys	Glu	Leu	Arg	Glu	Arg		
	515					520						525					
Val	Val	Lys	Glu	Leu	Asp	Lys	Val	Gly	Asn	Asn	Ala	Ser	Ala	Val	Ala		
	530					535					540						
Glu	Val	Asp	Phe	Ile	Gln	Ser	Leu	Ala	Gln	Ile	Ala	Tyr	Glu	Lys	Asp		
545					550					555					560		
Trp	Ala	Lys	Pro	Gln	Ile	His	Glu	Gly	Tyr	Glu	Leu	Ile	Ile	Glu	Glu		
				565					570					575			
Gly	Arg	His	Pro	Val	Ile	Glu	Glu	Phe	Val	Glu	Asn	Tyr	Val	Pro	Asn		
			580					585					590				
Asp	Thr	Lys	Leu	Asp	Arg	Asp	Ser	Phe	Ile	His	Val	Ile	Thr	Gly	Pro		
		595					600					605					
Asn	Met	Ala	Gly	Lys	Ser	Ser	Tyr	Ile	Arg	Gln	Val	Gly	Val	Leu	Thr		
	610					615					620						
Leu	Leu	Ser	His	Ile	Gly	Ser	Phe	Ile	Pro	Ala	Arg	Arg	Ala	Lys	Ile		
625					630					635					640		
Pro	Val	Val	Asp	Ala	Leu	Phe	Thr	Arg	Ile	Gly	Ser	Gly	Asp	Val	Leu		
				645					650					655			
Ala	Leu	Gly	Val	Ser	Thr	Phe	Met	Asn	Glu	Met	Leu	Glu	Val	Ser	Asn		
			660					665					670				
Ile	Leu	Asn	Asn	Ala	Thr	Glu	Lys	Ser	Leu	Val	Ile	Leu	Asp	Glu	Val		
	675						680					685					
Gly	Arg	Gly	Thr	Ser	Thr	Tyr	Asp	Gly	Ile	Ala	Ile	Ser	Lys	Ala	Ile		
	690					695				700							
Val	Lys	Tyr	Ile	Ser	Glu	Lys	Leu	Lys	Ala	Lys	Thr	Leu	Leu	Ala	Thr		
705					710					715					720		
His	Phe	Leu	Glu	Ile	Thr	Glu	Leu	Glu	Gly	Lys	Ile	Glu	Gly	Val	Lys		
			725						730					735			
Asn	Tyr	His	Met	Glu	Val	Glu	Lys	Thr	Pro	Glu	Gly	Ile	Arg	Phe	Leu		
			740					745					750				
Tyr	Ile	Leu	Lys	Glu	Gly	Lys	Ala	Glu	Gly	Ser	Phe	Gly	Ile	Glu	Val		
		755					760					765					
Ala	Lys	Leu	Ala	Gly	Leu	Pro	Glu	Glu	Val	Val	Glu	Glu	Ala	Arg	Lys		
	770					775					780						
Ile	Leu	Arg	Glu	Leu	Glu	Glu	Lys	Glu	Asn	Lys	Lys	Glu	Asp	Ile	Val		
785					790					795					800		
Pro	Leu	Leu	Glu	Glu	Thr	Phe	Lys	Lys	Ser	Glu	Glu	Ala	Gln	Arg	Leu		
			805						810					815			
Glu	Glu	Tyr	Glu	Glu	Ile	Ile	Lys	Lys	Ile	Glu	Glu	Ile	Asp	Ile	Gly		
			820					825					830				
Asn	Thr	Thr	Pro	Leu	Gln	Ala	Leu	Leu	Ile	Leu	Ala	Glu	Leu	Lys	Lys		
		835					840					845					

-continued

Lys Cys Ser Phe Ser Lys Lys Glu Ser Gly Ala
 850 855

<210> SEQ ID NO 5
 <211> LENGTH: 855
 <212> TYPE: PRT
 <213> ORGANISM: Aquifex pyrophilus

<400> SEQUENCE: 5

Met Gly Lys Glu Lys Lys Glu Leu Thr Pro Met Leu Ala Gln Tyr His
 1 5 10 15

Gln Phe Lys Ser Met Tyr Pro Asp Cys Leu Leu Leu Phe Arg Leu Gly
 20 25 30

Asp Phe Tyr Glu Leu Phe Tyr Glu Asp Ala Val Val Gly Ser Lys Glu
 35 40 45

Leu Gly Leu Val Leu Thr Ser Arg Pro Ala Gly Lys Gly Arg Glu Arg
 50 55 60

Ile Pro Met Cys Gly Val Pro Tyr His Ser Ala Asn Asn Tyr Ile Ala
 65 70 75 80

Lys Leu Val Asn Lys Gly Tyr Lys Val Ala Ile Cys Glu Gln Val Glu
 85 90 95

Asp Pro Ser Lys Ala Lys Gly Ile Val Lys Arg Asp Val Ile Arg Val
 100 105 110

Ile Thr Pro Gly Thr Phe Phe Glu Arg Glu Thr Gly Gly Leu Cys Ser
 115 120 125

Leu Tyr Arg Lys Gly Lys Ser Tyr Leu Val Ser Tyr Leu Asn Leu Ser
 130 135 140

Val Gly Glu Phe Ile Gly Ala Lys Val Lys Glu Glu Glu Leu Ile Asp
 145 150 155 160

Phe Leu Ser Lys Phe Asn Ile Arg Glu Val Leu Val Lys Lys Gly Glu
 165 170 175

Lys Leu Pro Glu Lys Leu Glu Lys Val Leu Lys Leu His Ile Thr Glu
 180 185 190

Leu Glu Glu Glu Phe Phe Glu Glu Gly Lys Glu Glu Leu Leu Lys Asp
 195 200 205

Tyr Gly Val Pro Ser Ile Lys Ala Phe Gly Phe Gln Asp Glu Asp Leu
 210 215 220

Ser Leu Ser Leu Gly Ala Val Tyr Arg Tyr Ala Lys Ala Thr Gln Lys
 225 230 235 240

Ser Phe Thr Pro Leu Ile Pro Lys Pro Lys Pro Tyr Val Asp Glu Gly
 245 250 255

Tyr Val Lys Leu Asp Leu Lys Ala Val Lys Gly Leu Glu Ile Thr Glu
 260 265 270

Ser Ile Glu Gly Arg Lys Asp Leu Ser Leu Phe Lys Val Val Asp Arg
 275 280 285

Thr Leu Thr Gly Met Gly Arg Arg Arg Leu Arg Phe Arg Leu Leu Asn
 290 295 300

Pro Phe Arg Ser Ile Glu Arg Ile Arg Lys Val Gln Glu Ala Val Glu
 305 310 315 320

Glu Leu Ile Asn Lys Arg Glu Val Leu Asn Glu Ile Arg Lys Thr Leu
 325 330 335

Glu Gly Met Ser Asp Leu Glu Arg Leu Val Ser Arg Ile Ser Ser Asn
 340 345 350

Met Ala Ser Pro Arg Glu Leu Ile His Leu Lys Asn Ser Leu Arg Lys

-continued

355	360	365
Ala Glu Glu Leu Arg Lys Ile Leu Ser Leu Leu Asp Ser Glu Ile Phe 370 375 380		
Lys Glu Ile Glu Gly Ser Leu Leu Asn Leu Asn Lys Val Ala Asp Leu 385 390 395 400		
Ile Asp Lys Thr Leu Val Asp Asp Pro Pro Leu His Val Lys Glu Gly 405 410 415		
Gly Leu Ile Lys Pro Gly Val Asn Ala Tyr Leu Asp Glu Leu Arg Phe 420 425 430		
Ile Arg Glu Asn Ala Glu Lys Leu Leu Lys Glu Tyr Glu Lys Lys Leu 435 440 445		
Lys Lys Glu Thr Gly Ile Gln Ser Leu Lys Ile Gly Tyr Asn Lys Val 450 455 460		
Met Gly Tyr Tyr Ile Glu Val Thr Lys Ala Asn Val Lys Tyr Val Pro 465 470 475 480		
Glu His Phe Arg Arg Gln Thr Leu Ser Asn Ala Glu Arg Tyr Thr 485 490 495		
Thr Glu Glu Leu Gln Arg Leu Glu Glu Lys Ile Leu Ser Ala Gln Thr 500 505 510		
Arg Ile Asn Glu Leu Glu Tyr Glu Leu Tyr Arg Glu Leu Arg Glu Glu 515 520 525		
Val Val Lys Glu Leu Asp Lys Val Gly Asn Asn Ala Thr Leu Ile Gly 530 535 540		
Glu Val Asp Tyr Ile Gln Ser Leu Ala Trp Leu Ala Leu Glu Lys Gly 545 550 555 560		
Trp Val Lys Pro Glu Val His Glu Gly Tyr Glu Leu Ile Ile Glu Glu 565 570 575		
Gly Lys His Pro Val Ile Glu Glu Phe Thr Lys Asn Tyr Val Pro Asn 580 585 590		
Asp Thr Lys Leu Thr Glu Glu Glu Phe Ile His Val Ile Thr Gly Pro 595 600 605		
Asn Met Ala Gly Lys Ser Ser Tyr Ile Arg Gln Val Gly Val Leu Thr 610 615 620		
Leu Leu Ala His Thr Gly Ser Phe Leu Pro Val Lys Ser Ala Arg Ile 625 630 635 640		
Pro Leu Val Asp Ala Ile Phe Thr Arg Ile Gly Ser Gly Asp Val Leu 645 650 655		
Ala Leu Gly Val Ser Thr Phe Met Asn Glu Met Leu Asp Val Ser Asn 660 665 670		
Ile Leu Asn Asn Ala Thr Lys Arg Ser Leu Ile Ile Leu Asp Glu Val 675 680 685		
Gly Arg Gly Thr Ser Thr Tyr Asp Gly Ile Ala Ile Ser Lys Ala Ile 690 695 700		
Val Lys Tyr Ile Ser Glu Lys Ile Gly Ala Lys Thr Leu Leu Ala Thr 705 710 715 720		
His Tyr Leu Glu Leu Thr Glu Leu Glu Arg Lys Val Lys Gly Val Lys 725 730 735		
Asn Tyr His Met Glu Val Glu Glu Thr Asp Glu Gly Ile Arg Phe Leu 740 745 750		
Tyr Ile Leu Lys Glu Gly Arg Ala Lys Gly Ser Phe Gly Ile Asp Val 755 760 765		
Ala Lys Leu Ala Gly Leu Pro Glu Glu Val Val Arg Glu Ala Lys Lys 770 775 780		

-continued

```

Ile Leu Lys Glu Leu Glu Gly Glu Lys Gly Lys Gln Glu Val Leu Pro
785              790              795              800

Phe Leu Glu Glu Thr Tyr Lys Lys Ser Val Asp Glu Glu Lys Leu Asn
      805              810              815

Phe Tyr Glu Glu Ile Ile Lys Glu Ile Glu Glu Ile Asp Ile Gly Asn
      820              825              830

Thr Thr Pro Val Lys Ala Leu Leu Ile Leu Ala Glu Leu Lys Glu Arg
      835              840              845

Ile Lys Ser Phe Ile Lys Arg
      850              855

```

```

<210> SEQ ID NO 6
<211> LENGTH: 793
<212> TYPE: PRT
<213> ORGANISM: Thermotoga maritima

```

```

<400> SEQUENCE: 6

```

```

Met Lys Val Thr Pro Leu Met Glu Gln Tyr Leu Arg Ile Lys Glu Gln
1              5              10              15

Tyr Lys Asp Ser Ile Leu Leu Phe Arg Leu Gly Asp Phe Tyr Glu Ala
      20              25              30

Phe Phe Glu Asp Ala Lys Ile Val Ser Lys Val Leu Asn Ile Val Leu
      35              40              45

Thr Arg Arg Gln Asp Ala Pro Met Ala Gly Ile Pro Tyr His Ala Leu
      50              55              60

Asn Thr Tyr Leu Lys Lys Leu Val Glu Ala Gly Tyr Lys Val Ala Ile
65              70              75              80

Cys Asp Gln Met Glu Glu Pro Ser Lys Ser Lys Lys Leu Ile Arg Arg
      85              90              95

Glu Val Thr Arg Val Val Thr Pro Gly Ser Ile Val Glu Asp Glu Phe
      100             105             110

Leu Ser Glu Thr Asn Asn Tyr Met Ala Val Val Ser Glu Glu Lys Gly
      115             120             125

Arg Tyr Cys Thr Val Phe Cys Asp Val Ser Thr Gly Glu Val Leu Val
      130             135             140

His Glu Ser Ser Asp Glu Gln Glu Thr Leu Asp Leu Leu Lys Asn Tyr
145             150             155             160

Ser Ile Ser Gln Ile Ile Cys Pro Glu His Leu Lys Ser Ser Leu Lys
      165             170             175

Glu Arg Phe Pro Gly Val Tyr Thr Glu Thr Ile Ser Glu Trp Tyr Phe
      180             185             190

Ser Asp Leu Glu Glu Val Glu Lys Ala Tyr Asn Leu Lys Asp Ile His
      195             200             205

His Phe Glu Leu Ser Pro Leu Ala Leu Lys Ala Leu Ala Ala Leu Ile
      210             215             220

Lys Tyr Val Lys Tyr Thr Met Ile Gly Glu Asp Leu Asn Leu Lys Pro
225             230             235             240

Pro Leu Leu Ile Ser Gln Arg Asp Tyr Met Ile Leu Asp Ser Ala Thr
      245             250             255

Val Glu Asn Leu Ser Trp Ile Pro Gly Asp Arg Gly Lys Asn Leu Phe
      260             265             270

Asp Val Leu Asn Asn Thr Glu Thr Pro Met Gly Ala Arg Leu Gly Lys
      275             280             285

Lys Trp Ile Leu His Pro Leu Val Asp Arg Lys Gln Ile Glu Glu Arg

```

-continued

290	295	300
Leu Lys Ala Val Glu Arg 305	Leu Val Asn Asp Arg Val Ser Leu Glu Glu 310	315 320
Met Arg Asn Leu Leu Ser Asn Val Arg Asp Val Glu Arg Ile Val Ser 325	330	335
Arg Val Glu Tyr Asn Arg Ser Val Pro Arg Asp Leu Val Ala Leu Arg 340	345	350
Glu Thr Leu Glu Ile Ile Pro Lys Leu Asn Glu Val Leu Ser Thr Phe 355	360	365
Gly Val Phe Lys Lys Leu Ala Phe Pro Glu Gly Leu Val Asp Leu Leu 370	375	380
Arg Lys Ala Ile Glu Asp Asp Pro Val Gly Ser Pro Gly Glu Gly Lys 385	390	395 400
Val Ile Lys Arg Gly Phe Ser Ser Glu Leu Asp Glu Tyr Arg Asp Leu 405	410	415
Leu Glu His Ala Glu Glu Arg Leu Lys Glu Phe Glu Glu Lys Glu Arg 420	425	430
Glu Arg Thr Gly Ile Gln Lys Leu Arg Val Gly Tyr Asn Gln Val Phe 435	440	445
Gly Tyr Tyr Ile Glu Val Thr Lys Ala Asn Leu Asp Lys Ile Pro Asp 450	455	460
Asp Tyr Glu Arg Lys Gln Thr Leu Val Asn Ser Glu Arg Phe Ile Thr 465	470	475 480
Pro Glu Leu Lys Glu Phe Glu Thr Lys Ile Met Ala Ala Lys Glu Arg 485	490	495
Ile Glu Glu Leu Glu Lys Glu Leu Phe Thr Ser Val Cys Glu Glu Val 500	505	510
Lys Lys His Lys Glu Val Leu Leu Glu Ile Ser Glu Asp Leu Ala Lys 515	520	525
Ile Asp Ala Leu Ser Thr Leu Ala Tyr Asp Ala Ile Met Tyr Asn Tyr 530	535	540
Thr Lys Pro Val Phe Ser Glu Asp Arg Leu Glu Ile Lys Gly Gly Arg 545	550	555 560
His Pro Val Val Glu Arg Phe Thr Gln Asn Phe Val Glu Asn Asp Ile 565	570	575
Tyr Met Asp Asn Glu Lys Arg Phe Val Val Ile Thr Gly Pro Asn Met 580	585	590
Ser Gly Lys Ser Thr Phe Ile Arg Gln Val Gly Leu Ile Ser Leu Met 595	600	605
Ala Gln Ile Gly Ser Phe Val Pro Ala Gln Lys Ala Ile Leu Pro Val 610	615	620
Phe Asp Arg Ile Phe Thr Arg Met Gly Ala Arg Asp Asp Leu Ala Gly 625	630	635 640
Gly Arg Ser Thr Phe Leu Val Glu Met Asn Glu Met Ala Leu Ile Leu 645	650	655
Leu Lys Ser Thr Asn Lys Ser Leu Val Leu Leu Asp Glu Val Gly Arg 660	665	670
Gly Thr Ser Thr Gln Asp Gly Val Ser Ile Ala Trp Ala Ile Ser Glu 675	680	685
Glu Leu Ile Lys Arg Gly Cys Lys Val Leu Phe Ala Thr His Phe Thr 690	695	700
Glu Leu Thr Glu Leu Glu Lys His Phe Pro Gln Val Gln Asn Lys Thr 705	710	715 720

-continued

Ile	Leu	Val	Lys	Glu	Glu	Gly	Lys	Asn	Val	Ile	Phe	Thr	His	Lys	Val
			725						730					735	
Val	Asp	Gly	Val	Ala	Asp	Arg	Ser	Tyr	Gly	Ile	Glu	Val	Ala	Lys	Ile
			740					745					750		
Ala	Gly	Ile	Pro	Asp	Arg	Val	Ile	Asn	Arg	Ala	Tyr	Glu	Ile	Leu	Glu
			755				760					765			
Arg	Asn	Phe	Lys	Asn	Asn	Thr	Lys	Lys	Asn	Gly	Lys	Ser	Asn	Arg	Phe
			770			775					780				
Ser	Gln	Gln	Ile	Pro	Leu	Phe	Pro	Val							
			785			790									

The invention claimed is:

1. A method for performing a multiplex single primer extension reaction, comprising:

a. extending at least 1,000 different primers in a single primer extension reaction using a plurality of target nucleic acids as templates in the presence of (1) a MutS protein, and (2) Mg^{2+} at a minimum concentration of 6 mM to generate extension products, wherein the primers each differ from the others.

2. The method of claim 1, wherein at least 4,000 additional different primers are extended in the single primer extension reaction.

3. The method of claim 1, wherein at least 9,000 additional different primers are extended in the single primer extension reaction.

4. The method of claim 1, wherein the single primer extension reaction is performed in the absence of a primer capable of specifically annealing to a portion of the extension products.

5. The method of claim 1, wherein most of the target nucleic acids comprise a common sequence that is at least 10 nucleotides in length, and wherein the reaction is performed in the presence of a boosting primer that comprises the common sequence or a portion thereof that is at least 10 nucleotides in length.

6. The method of claim 5, wherein each of the primers other than the boosting primer is present in the single primer extension reaction at a concentration of 1 nM to 10 nM.

7. The method of claim 1, wherein each of the different primers is present in the single primer extension reaction at a concentration of no more than 20 nM.

8. The method of claim 1, wherein Mg^{2+} is present in the single primer extension reaction at a concentration of 6 mM to 20 mM.

9. The method of claim 1, wherein the plurality of target nucleic acids are genomic DNA or amplification products of genomic DNA.

10. The method of claim 1, wherein the plurality of target nucleic acids are cDNA or amplification products of cDNA.

11. The method of claim 1, wherein the plurality of target nucleic acids are bisulfite treated genomic DNA or its amplified products.

12. The method of claim 1, wherein the MutS protein is a *Thermus aquaticus* MutS protein.

13. The method of claim 1, wherein the primer annealing and extension specificity of the single primer extension reaction is at least 80%.

14. The method of claim 1, wherein the average primer annealing and extension efficiency of the single primer extension reaction is at least 10%.

15. The method of claim 1, further comprising:

b. optionally amplifying the extension product(s) to generate amplified extension product(s), and

c. sequencing the extension product(s) or the amplified extension products.

16. The method of claim 1, wherein the target nucleic acids are partially double-stranded, and comprise a double-stranded target nucleic acid sequence and a single-stranded common sequence at the 5' terminus of each strand of the target nucleic acid sequence, wherein the single-stranded common sequence is at least 10 nucleotides in length.

17. The method of claim 16, comprising: prior to step a.,

(i) fragmenting DNA molecules to generate double-stranded DNA fragments, and

(ii) ligating the single-stranded common sequence to the 5' terminus of each strand of the double-stranded DNA fragment to provide partially double-stranded target nucleic acids.

18. The method of claim 1, wherein each of the different primers comprises at its 5' terminus a universal sequence that is incapable of specifically annealing to a region of a target nucleic acid.

19. A method for sequencing nucleic acids, comprising:

a. providing a plurality of partially double-stranded target nucleic acids each of which comprises a double-stranded target nucleic acid sequence, and a single-stranded common sequence at the 5' terminus of each strand of the target nucleic acid sequence, wherein the single-stranded common sequence is the same among the plurality of the target nucleic acids and is at least 10 nucleotides in length;

b. providing a plurality of primers capable of specifically annealing to portions of the target nucleic acid sequences;

c. extending the primers using the target nucleic acid sequences as templates in the presence of (1) a MutS protein and (2) Mg^{2+} at a minimum concentration of 6 mM to obtain extension products; and

d. sequencing the extension products or their amplification products.

20. The method of claim 19, wherein extending the primers in step c. is performed in the presence of a boosting primer comprising the common sequence or a portion thereof that is at least 10 nucleotides in length.

21. The method of claim 19, wherein step a. comprises:

a1. fragmenting double-stranded DNA molecules to generate double-stranded DNA fragments, and

a2. ligating the single-stranded common sequence to the 5' terminus of each strand of the double-stranded DNA fragments to provide the plurality of partially double-stranded nucleic acids.

22. The method of claim **19**, comprising prior to step d., a further step of amplifying the extension products of step c.

23. A kit for performing a primer extension reaction, comprising:

- (1) a pool of at least 1,000 different primers, wherein the 5 primers each differ from the others,
- (2) a buffer concentrate comprising a MutS protein, dNTPs, and Mg^{2+} , wherein the concentration of Mg^{2+} in a highly multiplex single primer extension reaction mixture generated from diluting the buffer concentrate 10 is at least 6 mM, and
- (3) a DNA polymerase.

24. The kit of claim **23**, wherein the kit comprises a pool of at least 10,000 different primers.

25. The kit of claim **23**, wherein the MutS protein is a 15 *Thermus aquaticus* MutS protein.

* * * * *