

ClearCorrect Support > FAQ > Aligner FAQs

Aligner Material - FAQs

Updated 2 years ago

Table of contents

What are ClearCorrect aligners made out of?

[Are ClearCorrect aligners BPA-free?](#)

[Can my patient drink hot beverages like coffee or tea while wearing aligners?](#)

[Is the Zendura used to make the aligners the same as the Zendura sold in the store?](#)

[What machines can I use to make my own Zendura aligner?](#)

[What plastic do you recommend for making replacement aligners?](#)

[Why are the distal edges of the aligner missing sometimes?](#)

For frequently asked questions on aligner material, refer to the following:

- [What are ClearCorrect aligners made out of?](#)
- [Are ClearCorrect aligners BPA-free?](#)
- [Can my patient drink hot beverages like coffee or tea while wearing aligners?](#)
- [Is the Zendura used to make the aligners the same as the Zendura sold in the online store?](#)
- [What machines can I use to make my own Zendura aligner?](#)
- [What plastic do you recommend for making replacement aligners?](#)
- [Why are the distal edges of the aligner missing sometimes?](#)

ALIGN - EXHIBIT 2022
ClearCorrect Operating, LLC. v. Align Technology, Inc.
IPR2025-00821

For further information about our aligner material and safety information, please [contact your local Customer Care](#) representative.

What are ClearCorrect aligners made out of?

Our aligners are made from [Zendura](#), a custom formulation of medical-grade polyurethane which has been thoroughly tested for biocompatibility. This advanced plastic was created exclusively for ClearCorrect by [Bay Materials](#).

[Back to Top](#)

Are ClearCorrect aligners BPA-free?

Yes! Our aligners are FDA approved and free of bisphenol-A ("BPA") and phthalates.

[Back to Top](#)

Can my patient drink hot beverages like coffee or tea while wearing aligners?

We don't recommend it. Hot liquids can distort the material. Coffee and tea are also likely to stain the aligners, and sugar can get trapped between the teeth and the aligners. It's best to remove the aligners before eating or drinking anything other than water.

[Back to Top](#)

Is the Zendura used to make the aligners the same as the Zendura sold in the store?

Yes! The same Zendura plastic used to make your patients' aligners is [available on our online store](#) at a providers-only price of \$60 for 20 sheets. This way, you won't have to compromise quality for the convenience of making your own aligners in-house.

[Back to Top](#)

What machines can I use to make my own Zendura aligner?

[Bay Materials recommends](#) the use of a pressure-forming machine (such as a Biostar, MiniStar, or Druformat) when thermoforming Zendura plastic. However, you may be able to achieve good results using a vacuum-forming machine by following [these instructions](#).

[Back to Top](#)

What plastic do you recommend for making replacement aligners?

ClearCorrect aligners are made from **Zendura®** plastic which has been specially formulated for ClearCorrect by [Bay Materials](#). It is a superior plastic for clear aligners: it leads the class in stress retention, crack & impact resistance, clarity, and stain resistance.

Zendura aligner plastic is [available for purchase](#) at a special discounted price for ClearCorrect providers in the online store.

Consult your thermoforming machine manual for proper instructions, heating time, cooling time, etc.

[Back to Top](#)

Why are the distal edges of the aligner missing sometimes?

Sometimes we receive impressions that don't have enough detail to accurately model the distal edges of the posterior teeth. Distortion is much more prevalent in this area, because it can be difficult to make sure that the impression material completely covers the teeth in the back of the mouth.

Of course, we always prefer to receive complete, accurate impressions. But we don't want to inconvenience our providers unnecessarily either.

In some cases, we'll make an exception and process the case even though the distal surfaces of the posterior teeth are incomplete in our model. The aligners still have plenty of surface area to grip the teeth. We just trim off the potentially inaccurate area so that the case can progress without delay.

If you want to make sure that your patient's aligners fully cover the distal surfaces of all the teeth, just double-check your impressions to make sure that they're not distorted in that area, and you should be fine.

[Back to Top](#)

Articles in this section

Aligner Material - FAQs

ClearCorrect® Aligners FAQs

Case & Pricing Packages - FAQs

Was this article helpful?

Yes

No

25 out of 27 found this helpful

Comments

0 comments


Please [sign in](#) to leave a comment.

 Search...



Have more questions?

Submit a request



Share it, if you like it.



ClearCorrect	+
Products and Solutions	+
Services	+
Highlights	+
Straumann Group	+
Straumann Group Focus	+

