

< Edit Treatment



12:33 AM

Guest User



Signature Hydrafacial Exfoliation

Instructions

Adjust the serum flow valve on your handpiece as needed. Adjust vacuum setting as needed. Use a dab and spread motion prior to the first pass to moisten skin. Use long, fast, 60% overlapping strokes. Make two passes for complete exfoliation. Hold skin taut and apply counter pressure.

Treatment Pattern

○ Dab ◄ Long, Fast, Overlapping Strokes



Step 1



Step 2



Step 3



Step 4

Vacuum



25



START

Device



Handpiece

Tip



Blue

Serum



Activ-4



ACTIV-4

BETA HD CLEAR

HYDRATE

RINSEAWAY



Aesthetic Management Partners,
LLC et al v. HydraFacial
IPR2025-01217
HydraFacial Ex. 2128