
From: Nathan Moskowitz <nmosk@hotmail.com>
Sent: Monday, August 8, 2011 11:53 PM CDT
To: eric sugalski <eric@bostondevice.com>
CC: andy moskowitz <andymoskow@gmail.com>
Subject: re:screwdrivers

Eric,

I just received the 2 screwdrivers from Mcmasters today. They both fit into the prototype and easily open and close the apparatus. The mechanism works beautifully.

Great job.

Thanks

Nathan

From: Eric Sugalski <eric@bostondevice.com>
Sent: Tuesday, August 9, 2011 7:48 AM CDT
To: Nathan Moskowitz <nmosk@hotmail.com>
CC: andy moskowitz <andymoskow@gmail.com>; David Schoon <dschoon@bostondevice.com>
Subject: Re: screwdrivers

Great to hear. Thank you for the feedback, Nathan. Sorry we didn't get you the correct screwdrivers earlier!

Eric

--

Eric Sugalski
Boston Device Development, Inc
office - (617) 279-0057 x703
mobile - (617) 852-6622
www.bostondevice.com

On Tue, Aug 9, 2011 at 12:53 AM, Nathan Moskowitz <nmosk@hotmail.com> wrote:

Eric,

I just received the 2 screwdrivers from Mcmasters today. They both fit into the prototype and easily open and close the apparatus. The mechanism works beautifly.

Great job.

Thanks

Nathan

From: Eric Sugalski <eric@bostondevice.com>
Sent: Friday, August 5, 2011 9:47 AM CDT
To: Andy Moskowitz <andymoskow@gmail.com>
Subject: Re: checking in

Andy - Apologies for this. We are sending you some new drivers. You should have them early next week.
Eric

Sent from my iPhone

On Aug 5, 2011, at 8:11 AM, Andy Moskowitz <andymoskow@gmail.com> wrote:

Dear Eric,

My father said that he tries the parts we got

[7996A83](#) and

[6286A12](#) and it did not fit is there anything else we can try?

Andy

On Sat, Jul 23, 2011 at 10:49 PM, Eric Sugalski <eric@bostondevice.com> wrote:

Hi Andy & Nathan -

Just checking in regarding the prototype. Have you had a chance to evaluate it?

Let me know your thoughts when you get a chance.

Thanks,

Eric

--

Eric Sugalski
Boston Device Development, Inc
office - [\(617\) 279-0057 x703](tel:(617)279-0057)
mobile - [\(617\) 852-6622](tel:(617)852-6622)
www.bostondevice.com

From: Nathan Moskowitz <nmosk@hotmail.com>
Sent: Monday, July 25, 2011 10:34 PM CDT
To: eric sugalski <eric@bostondevice.com>
Subject: RE: checking in

Thanks Eric

From: eric@bostondevice.com
Date: Mon, 25 Jul 2011 22:23:01 -0400
Subject: Re: checking in
To: nmosk@hotmail.com
CC: andymoskow@hotmail.com

Hi Nathan -

Thanks for getting back to me. There are a wide variety of drivers for the M2 allen used in that device. Following are a few part numbers from McMaster.com that should work:

[7996A83](#)
[6286A12](#)

Let me know if you think either of these will do the trick.

Eric

--

Eric Sugalski
Boston Device Development, Inc
office - (617) 279-0057 x703
mobile - (617) 852-6622
www.bostondevice.com

On Sun, Jul 24, 2011 at 12:00 PM, Nathan Moskowitz <nmosk@hotmail.com> wrote:

Hi Eric,

It looks very nice.

Unfortunately I can't evaluate its functional utility of expansion because I don't have any instrument to fit into the very small screw insert.

When functionally placing this in the spine, a surgeon can not use an allen wrench because you'll be bumping into nerves and paralyzing patients. What is needed is a straight screw driver instrument with a distal protrusion that can mate with and interact with the insert which will allow expansion and collapse of the cage.

Please inform us where we can purchase such an instrument with the precisediameter needed to couple with your insert.

Thanks,
Nathan

From: eric@bostondevice.com
Date: Sat, 23 Jul 2011 21:49:24 -0400
Subject: checking in
To: andymoskow@gmail.com; nmosk@hotmail.com

Hi Andy & Nathan -
Just checking in regarding the prototype. Have you had a chance to evaluate it? Let me know your thoughts when you get a chance.
Thanks,
Eric

--
Eric Sugalski
Boston Device Development, Inc
office - (617) 279-0057 x703
mobile - (617) 852-6622
www.bostondevice.com