


NXT Hg5[®] Series by Solmetex[™]



The NXT Hg5 is the next generation in amalgam separation and the gold standard. Our products meet or exceed all state and federal requirements for the safe removal of mercury in the dental office.

Our signature product – the NXT Hg5 comes in several sizes to fit the needs of your practice.



NXT Hg5


NXT Hg5 mini

Most Popular

Small Office

[Learn more](#)

[Learn more](#)



NXT Hg5 High Volume

Large Office

[Learn more](#)

Recycling Kit

Easy Disposal Solution

[Learn more](#)

The proven answer to your amalgam separation concerns.

Solmetex provides the complete solution:

- ISO 11143:2008 certified by NSF International.
- Professionally engineered, robust, simple design.
- Easy to operate and maintenance free.
- Flexible design for hard-to-fit applications.
- Functional for wet or dry vacuum systems.





[Click Here](#) to read our Product Evaluation!

NXT Hg5

- Accommodates up to 10 operatories.
- Dimensions 30"h x 11"w x 8"d.
- 1 1/2" inlet and outlet piping.
- Fernco couplers for easy install.
- Floor/wall mount brackets provided.

Description and Part Number

- NXT Hg5 Amalgam Separator (NXT-HG5-001)
- NXT Hg5 Collection Container with Recycle Kit (NXT-HG5-002CR)

Support Data

- NXT Hg5 Series Brochure
- NXT Hg5 Installation & Operation Manual
- NXT Hg5 Installation Video





Order Now



NXT Hg5 mini

- Accommodates up to 4 operatories.
- Dimensions 25"h x 11"w x 8"d.
- 3/4" inlet and outlet piping.
- Fernco couplers for easy install.
- Floor/wall mount brackets provided.
- Designed specifically to fit under a counter.

Description and Part Number

- NXT Hg5 mini Amalgam Separator (NXT-HG5-MINI)
- NXT Hg5 Collection Container with Recycle Kit (NXT-HG5-002CR)

Support Data

- NXT Hg5 Series Brochure
- NXT Hg5 mini Installation & Operation Manual
- NXT Hg5 Installation Video



Order Now

NXT Hg5 High Volume

- Accommodates up to 20 operators.
- Dimensions 29 3/4" h x 17 5/8" w x 22 1/2" d.
- 3" inlet and outlet piping.
- Wall mount brackets provided.

Description and Part Number

- NXT Hg5 High Volume Amalgam Separator (NXT-Hg5-HV)
- Collection Container with Recycle Kit (NXT-Hg5-002CR)

Support Data

- NXT Hg5 High Volume Series Brochure
- NXT Hg5 High Volume Installation & Operation Manual



Order Now

NXT Hg5 Series of Collection Containers with Recycle Kits

Once you are set up with one of our award winning NXT Hg5 series of amalgam separators that fits the size of your practice, in time, you will need to swap out the collection container located at the bottom of the amalgam separator for recycling. The frequency with which you change the collection container depends upon the size and volume of your practice, the container will need to be replaced once the sediment reaches the full line or once a year (whichever comes first). The good news is, we do all of the hard work and dispose of the mercury under federal guidelines for you. All you have to do is disconnect the old collection container and ship it in the provided recycling kit box to our recycling facility. Simply insert the new collection container, and you're ready to go. Keeping harmful mercury out of the environment has never been easier!

Description

- NXT Hg5, Hg5 [Legacy] or RAMVAC Utility Collection Container
- Solmetex box for shipping & recycling
- Shipping label
- Recycling documentation

Support Data

- NXT Hg5 Collection Container Replacement & Recycle Procedure
- NXT Shipping Code & Label Instructions
- NXT Hg5 Collection Container Change-Out Video
- Hg5 Collection Container Replacement & Recycle Procedure
- Hg5 Shipping Code & Label Instructions





Hg5 Collection Container

Order Now



NXT Hg5 Collection Container

Order Now



RAMVAC Hg5 Collection Container

Order Now



Solmetex™, LLC

English

4 Mount Royal Ave, Ste 200

Marlborough, MA 01752-1983

(800) 216-5505 Main

(508) 393-1795 Fax

Hours: 8:30am-8:00pm EST

Join Our Newsletter

* Email



By submitting this form, you are consenting to receive marketing emails from: Solmetex, 50 Bearfoot Rd, Northborough, MA, 01532, US, <https://solmetex.com/>. You can revoke your consent to receive emails at any time by using the SafeUnsubscribe® link, found at the bottom of every email. Emails are serviced by Constant Contact.

Subscribe

Partnerships





[Careers](#) • [Press Room](#) • [Privacy Policy](#)

Copyright © 2023 Solmetex

