


**MiniMed Paradigm®
REAL-Time System**
Guardian® REAL-Time

Benefits

► Features

Components

CareLink™ Personal Software

Skins

Pediatric Model

Tutorial

Insurance Coverage

**CareLink™ Therapy
Management Software**
Infusion Sets
CGMS® System Gold™
Other Insulin Pumps
How to Learn More
**Request Insulin
Pump Info**
**Request Guardian
REAL-Time Info**
**Patient Information
Events**
Find a Physician
**Important Safety
Information**
E-mail This Page
Print This Page
[Home](#) ► [Product Information](#) ► [Guardian REAL-Time](#) ► [Features](#)

Features that fit your diabetes management lifestyle

Type Size **+** **-**

The next level...



Want to go to the next level with the Guardian REAL-Time System?

[Try the Guardian Tutorial>](#)

As the world leader in continuous glucose monitoring, Medtronic Diabetes is committed to providing our customers with unique value that goes beyond delivering state-of-the-art technology. Your Guardian® REAL-Time Continuous Glucose Monitoring System includes:

- REAL-Time display for added protection, 24 hours a day.
- MiniLink™ REAL-Time Transmitter designed for exceptional comfort, with an ultra-sleek, thin profile.
- CareLink™ Personal software to help you manage your diabetes.
- Tools and support when you need them 24/7/365 from the most personalized, in-depth education and support program in the industry — by your side every step of the way.
- Cool Medtronic "Skins" for a personalized look.

REAL-Time display

The Guardian REAL-Time System display gives you added protection and confidence, 24 hours a day.



With the Guardian REAL-Time System, an amazing amount of information is right there on your monitor.

REAL-Time Alerts



Three different kinds of alerts give you time to avoid or minimize oncoming highs and lows.

REAL-Time Continuous Glucose Readings



When you can read your current glucose level any time of day or night, you won't be surprised by a sudden low or high.

REAL-Time Trend Graphs



Seeing a graph of your glucose levels over the last 3, 6, 12, or 24 hours lets you identify patterns and fine-tune your diabetes management accordingly.

REAL-Time Trend Arrows



When your glucose levels start changing in a hurry, these arrows show you the direction and rate of change, which helps you anticipate and react to the changes.

MiniLink REAL-Time Transmitter

The MiniLink REAL-Time Transmitter is exceptionally comfortable because it is small, lightweight, ergonomically shaped, and designed to minimize skin irritation. An ultra-sleek, thin profile allows the MiniLink REAL-Time Transmitter to lay flat without a bulge.

- Small and thin profile (approximately 1/3 of an inch thick) — easy to hide under clothing.
- Wireless — provides flexibility for placement in different locations on the body.
- Waterproof — can be worn while swimming*, bathing, or showering without worries

* Can be worn to a depth of 8 feet for 30 minutes

DEXCOM EXHIBIT 1056
Dexcom v. Abbott Diabetes
IPR2023-01409



Ultra-sleek design requires very little skin "real estate."


By your side every step of the way

When you start living with your Guardian REAL-Time System, you'll have the support of Medtronic Diabetes at your fingertips. A live, 24-Hour Helpline is available 365 days a year, including weekends and holidays. Highly-trained technicians answer your calls personally — no paging or forwarding. You can count on us to provide you with the industry's most personalized, in-depth education and support program.



Better control with CareLink Personal software

Web-based tools enable you to capture continuous glucose monitoring data and logbook entries online, and download easy-to-read charts, graphs, and tables. You and your healthcare provider can use these reports to discover patterns, problems, and insights not possible with fingerstick testing and logbooks. These insights can allow your healthcare provider to make therapy adjustments that can lead to better management of your diabetes. The reports can also help you identify opportunities for lifestyle changes that can improve control of your glucose levels.



Express yourself with Medtronic "Skins"

You can enhance the look of your Guardian REAL-Time System with "skins" — replaceable vinyl covers for your monitor that are a fun and exciting way to express your own style. Each high quality vinyl skin has an adhesive side, allowing you to easily apply the skin to your Guardian REAL-Time System monitor.

Go Ahead. [Express Yourself](#) today



[See the main components of the Guardian REAL-Time System](#)


[Store](#) | [CareLink](#) | [Pump School Online](#) | [Local Events](#) | [Contact Us](#)
[Home](#) | [Insulin Pump Therapy](#) | [Product Information](#) | [Help & Support](#) | [Healthcare Professionals](#) | [About Us](#)

MiniMed Paradigm® REAL-Time System

Guardian® REAL-Time

[Benefits](#)
[Features](#)
[Components](#)
[CareLink™ Personal Software](#)
[Skins](#)
[Pediatric Model](#)
[Tutorial](#)
[Insurance Coverage](#)

CareLink™ Therapy Management Software

Infusion Sets

CGMS® System Gold™

Other Insulin Pumps

[How to Learn More](#)
[Request Insulin Pump Info](#)
[Request Guardian REAL-Time Info](#)
[Patient Information Events](#)
[Find a Physician](#)
[Important Safety Information](#)
[E-mail This Page](#)
[Print This Page](#)
[Home](#) > [Product Information](#) > [Guardian REAL-Time](#) > [Components](#)

Components

Type Size [+](#) [-](#)

Want to learn more...



Want to learn more about the parts of the Guardian REAL-Time System?

[Try the Guardian Tutorial>](#)

Your Guardian® REAL-Time Continuous Glucose Monitoring System is compact, comfortable, and easy to wear and use. The small, comfortable, and discreet rechargeable transmitter connects directly to the glucose sensor and sends data wirelessly to the monitor, which can be up to 6 feet away. And the MiniLink™ REAL-Time Transmitter and glucose sensor are waterproof, allowing you to shower, bathe, and engage in water activities. ¹



The MiniLink REAL-Time Transmitter wirelessly transmits data collected by the glucose sensor to the monitor.

The Guardian REAL-Time System includes:

- Monitor
- MiniLink REAL-Time Transmitter and charger
- Glucose Sensor

Monitor

Your Guardian REAL-Time System monitor is small, discreet, and easy to wear — it can be attached to your belt, hidden in your pocket, or placed under your clothing.



MiniLink REAL-Time Transmitter and Charger

So small it makes a big difference — the MiniLink REAL-Time Transmitter is a key component of the Guardian REAL-Time System. Exceptionally small, lightweight, and ergonomically shaped, the transmitter can be worn while swimming, bathing or showering without worries. ¹ Its tiny rechargeable battery allows for 3 days of continuous use on a 20-minute recharge, and up to 14 days of continuous use on a full recharge.



The MiniLink REAL-Time Transmitter kit includes the transmitter, charger, and tester.

Glucose Sensor

The tiny glucose sensor is easily inserted, virtually painlessly, under the skin using an autoinsertion device that comes with the system.



¹ While the Guardian REAL-Time Continuous Glucose Monitoring System is water-resistant and designed to protect against splashes and occasional dunking, Medtronic Diabetes recommends that no continuous glucose monitoring system should be intentionally submerged in water.