


[Atraumatic Inserts](#)[Clamps & Clips](#)[Cosgrove Flex Clamps](#)[Occlusion Catheter](#)[Related Products](#)[Clot Management](#)[Carotid Shunts](#)[Biliary](#)[Product Catalog](#)[Request Information](#)[Reimbursement](#)[Products](#) / [Vascular Surgery](#) / [Atraumatic Occlusion](#) / [Occlusion Catheter](#)[Brochure](#)[Models](#)

Fogarty Occlusion Catheter

Occlusion catheters are available with a small balloon or a large balloon for use in a variety of procedures.

Small balloon occlusion catheters are effective in gaining temporary occlusion of blood flow in the following applications:

- Small diameter vessels such as those in the distal periphery.
- Operative occlusion of blood supply to specific organs such as the kidneys.
- Vascular-neuro surgery procedures, such as carotid-cavernous sinus fistulas.

Large balloon occlusion catheters are ideally suited for occlusion control in:

- Aortic or cardiac hemorrhage.
- Internal jugular vein hemorrhage near the base of the skull.
- Operation on the lower part of the abdominal aorta and renal arteries to eliminate the need for extensive dissection.
- Thrombectomy procedure utilizing the Fogarty graft thrombectomy catheter.

CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

[Contact Us](#) | [Legal Terms](#) | [Privacy Policy](#) | © 2015 Edwards Lifesciences Corporation. All rights reserved.